



The One That Works!
Disposable Air/Water Syringe Tips




Superior Spray
Instant Separation
Dry Air



Features

- » No special adapter required.
- » Interchangeable with metal tips.
- » Instant air/water separation.
- » No Leaks!
- » Superior air/water spray.
- » Beveled edges for patient safety.
- » Made from FDA approved plastic.
- » Guaranteed Dry Air!

Unique Quad Air Core Design



ATTENTION METAL TIP USERS

Are you aware that recent studies show metal tips are actually less sanitary than disposables?

That's because the very fine lumen, or air channels, actually act as an air trap, making the lumen susceptible to biofilm coating and collecting micro-debris from every patient. Even autoclaved sterilization can't clear these narrow channels of the debris that continue to build up and either clog your lines, or worse -- get blown into the mouth of your next patient!

Biofilm buildup also prevents condensation of steam on the walls of the lumen, hence making sterilization very difficult.



Rainbow Assorted



Lavender Purple



Marine Green



Polar White



Sunflower Yellow



Pacific Blue



Lilac Purple



Tangerine Orange



Orchid Purple



Mint Green



1-800-456-3660
stardentalsupply.com



Only Crystal Tip offers unique rewards, from gift certificates to Nintendo Wiis, and more. And when you make the switch to Crystal Tip, it is pain free. You receive free adapters (if required) and a great support system. And, only Crystal Tip offers a \$250 Working Guarantee. Please visit crystaltip.com for complete information, and instructions on how to convert your syringe from a proprietary adapter to its original manufactured condition so that you can use Crystal Tip properly.

MADE IN THE USA