



A Guide To Surgical Sutures






LOOK[®]



Absorbable

Material	Material Description	BSR Profile	Mass Absorption
Gut  • Surgical Gut Suture - Plain/Chromic * Plain Gut/Chromic Gut	Twisted • Plain & Chromic • Packaged Dry or in Wetting Solution Common uses: General Closure, OB/Gyn, Urology, Bowel Anastomosis	In vivo strength retention: Varies	Varies depending upon type of suture (plain vs. chromic) and site of placement
PolySyn FA[™] (glycolic acid)  • VICRYL RAPIDE [™] * Velosorb [™] Fast	Braided • Undyed (clear) • Coated Common uses: General Closure, OB/Gyn, Urology	In vivo strength retention: 48% at 7 days 0% at 14 days	Essentially complete between 49 and 63 days
PolySyn[™] (glycolic acid)  • VICRYL [™] * Polysorb [™] / Dexon [™]	Braided • Undyed (clear) • Dyed (violet) Monofilament • Undyed (clear) • Dyed (violet) • Coated & Uncoated Common uses: Skin Closure, OB/Gyn, Urology, Gastrointestinal	In vivo strength retention: 50% at 14 days 20% at 21 days	Essentially complete between 50 and 90 days

Non-Absorbable

Material	Material Description	Qualities
Silk  • PERMA-HAND [™] Silk * SofSilk [™]	Braided • Undyed (white) • Dyed (black) • Coated Common uses: General & Skin Closure, Gastrointestinal, Cardiovascular Surgery, Plastic Surgery, Ophthalmology	• Excellent handling and tying characteristics
Nylon  • ETHILON [™] * Monosof [™] /Dermalon [™]	Monofilament • Undyed (clear) • Dyed (blue) (black) • Packaged dry or in wetting solution Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery, Gastrointestinal	• Passes easily through tissue
Polypropylene  • PROLENE [™] * Surgipro [™] /Surgipro [™] II	Monofilament • Undyed (clear) • Dyed (blue) Common uses: General & Skin Closure, Neurosurgery, Plastic Surgery, Ophthalmology, Microsurgery, Cardiovascular Surgery	• Inert • Easy to handle • Passes easily through tissue
Polyviolene[™] (Uncoated) (poly (ethylene-terephthate))  • MERSILENE [™] * Surgidac [™]	Braided • Undyed (white) • Dyed (green) • Uncoated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security
Polyviolene[™] (Coated) (poly (ethylene-terephthate))  • ETHIBOND EXCEL [™] * Ti-Cron [™]	Braided • Undyed (white) • Dyed (green) • Coated Common uses: General & Skin Closure, Plastic Surgery, Ophthalmology, Cardiovascular Surgery, Neurological Procedures	• Smooth • Easy to handle • Good knot security

• Comparable Ethicon[™] Material * Comparable Syneture[™] Material

Cold Packs

(Graham)
Icepack 9100A
6.5" x 8.5", 24/case
855-0550 \$24.95



Grafco 010210
5" x 7", 48/case
855-0560 \$31.95



Coldstar 010407
5" x 5", 80/case
855-0570 \$45.95



(Allegiance)
Disposable Single-Use Packs; Cold Pack
activated by squeezing. Lasts 25 minutes.

Instant Cold Pack 10104
6 1/4" x 8", 16/case
855-0500 \$20.95



Soft-Pack 11445-020
6" x 6 1/2", 32/case
855-0510 \$26.95

Kwik-Kold 103B
4" x 5 3/4", 50/case
855-0520 \$32.95



(Nortech)
DentalPack®
Cold compress, good for up to 20 min-
utes, Reusable and Latex Free
4" x 5.75"
100/pack 855-0579 \$79.95



Periodontal Dressing

Periocare® (Pulpdent)
2 Paste System, perio Dressing, sets
resiliently hard. Won't chip or fall apart.

Kit 855-1015 \$48.95
3 oz tube paste, 3 oz tube gel & mixing
pad



Coe-Pak (GC)
Regular set is ready in 3 minutes, Hard
& Fast set is ready in 1 minute.

Standard pkg \$94.95
Regular Set 135001 855-0010
Fast/Hard Set 135301 855-0011



Coe-Pak Automix NDS (GC)
2-Double Barrel Cartridges (50ml
each), 12 mixing tips, instructions

855-0015 \$117.95



Noneugenol Surgical Dressing and Periodontal Pack
Leader of the Pack



COE-PAK™ AUTOMIX NDS
Periodontal Dressing

COE-PAK AUTOMIX NDS is a noneugenol surgical
dressing and periodontal pack that has no burning sensation,
no unpleasant taste, no discolorable odor, and offers proven
protection to surgical sites.

- Resilient Hardness
- Contains No Eugenol
- New Delivery System
- Two Year Shelf Life
- No Offensive Taste or Odor

The NDS version, New Delivery System, eliminates manual
mixing and provides a consistent mix for immediate
placement.



GC
GC AMERICA INC.

* Prices are subject to change without notice.



SURGICAL

Safety Scalpels

sterile

Uniblade™ (Applied Dental)
Disposable stainless steel scalpels
10 individually packs per box

#10	865-2200	\$10.95
#11	865-2201	\$7.95
#15	865-2204	\$7.95



Disposable Scalpels (Feather)

20/pkg, sterile \$19.95

#10	865-2010	#15C	865-2026
#11	865-2011	#20	865-2020
#15	865-2015	#21	865-2021



Feather Safeshield (Feather)
Disposable Sterile Scalpels

10/pkg


#10	865-2110	\$23.95
#11	865-2111	\$23.95
#15	865-2115	\$23.95
#15C	865-2118	\$26.95
#25	865-2125	\$23.95



Fax your order to: 310.537.1502

Surgical Blades

sterile



Stainless Steel
100/box

\$12.95

#10	865-0910	#15	865-1001	#22	865-0922
#11	865-0911	#15C	865-0918	#23	865-0923
#12	865-1000	#20	865-0920	#25	865-0925
#12B	865-0914	#21	865-0921		

Carbon Steel (Demetech)
100/box


\$11.95

#10	865-0110	#15	865-0115	#22	865-0122
#11	865-0111	#15C	865-0118	#23	865-0123
#12	865-0112	#20	865-0120	#25	865-0125
#12B	865-0114	#21	865-0121		

Stainless Steel (Demetech)
100/box


\$19.95

#10	865-0210	#15	865-0215	#22	865-0222
#11	865-0211	#15C	865-0218	#23	865-0223
#12	865-0212	#20	865-0220	#25	865-0225
#12B	865-0214	#21	865-0221		



Blade Handles

Blade Handle Regular (FCI)
 Each \$3.95
 Regular 865-8000
 With ruler 865-8001



(Miltex)
 Blade Handle # 3, 4-7R
 For blade sizes 10, 11, 12, 12b, 15 & 15C



865-8115 \$7.95

Blade Remover

Safety Blade Remover (Graham-Field)
 Individually packed, sterile
 10/pkg 875-7000 \$29.95



Blade Remover (Graham-Field)
 with container, holds between
 150-300 used blades
 Each 875-7010 \$5.95



Blade Remover (Integra/Miltex)
 Safely removes blades from handles,
 protecting hands, 4-100
 875-7020 \$17.95



Dry Socket Paste

Paste (Sultan)
 1 oz 857-4000 \$139.95



Surgical Pack

O-Pak® (Applied Dental)
 Sterile drape pack, with 14 essential
 disposable products including gowns
 and drapes. Latex free, sterilized by EO
 gas. 6 sealed packs per box



857-1500 \$68.95

Irrigating Syringes

412 Syringe
 12ml curved tips
 50/pkg

(Vista) 857-5000 \$22.95
 (Pac-Dent) 857-5010 \$19.95



Monoject Syringe #412
 (vMedtronic)

12ml, 50/pkg 857-5550 \$25.95



SURGICAL

Sutures

silk

Non-Absorbable, Silk, black braided (Look)



3/8 Circle Reverse Cutting, 12/pkg

784B	C-6	3-0	875-4700	\$27.95
781B		4-0	875-4701	\$26.95
754B		5-0	875-4702	\$41.95
786B	C-7	3-0	875-4703	\$25.95
770B	C-8	2-0	875-4704	\$25.95
787B	C-17	4-0	875-4705	\$25.95

778B	C-26	3-0	875-4750	\$25.95
785B	C-31	3-0	875-4751	\$25.95

(Sharp)
Non-Absorbable, Silk, black braided
3/8 Circle Reverse Cutting 12/pkg

784B	C-6	3-0	875-5001	\$11.95
------	-----	-----	----------	---------



1/2 Circle Reverse Cutting, 12/pkg

785B	C-31	3-0	875-5012	\$11.95
------	------	-----	----------	---------

(Applied Dental)
Unify Silk \$15.95
3/8 circle Reverse Cutting 18" long
12/pkg

Medium - 19mm needle

PS-2/FS-2	3-0	875-2000
	4-0	875-2001
	5-0	875-2002

Medium - 16mm needle

PS-3	4-0	875-2010
	5-0	875-2011



Non-Absorbable, Silk, black braided
(Covidien)



36/pkg	C-13	3-0	875-3300	\$59.95
12/pkg	C-13	4-0	875-3301	\$24.95
36/pkg	C-23	3-0	875-3302	\$59.95
36/pkg	C-1	5-0	875-3303	\$59.95

Sutures

plain gut

Absorbable, Plain Gut (Look)
3/8 Circle Reverse Cutting, 12/pkg

596B	C-3	5-0	875-4800	\$69.95
591B	C-6	4-0	875-4801	\$30.95
594B	C-31	3-0	875-4802	\$29.95



(Sharp)
Absorbable, Plain Gut \$19.95
3/8 Circle Reverse Cutting, 12/pkg

594B	C-31	3-0	875-5021	\$19.95
------	------	-----	----------	---------



(Applied Dental)
Premium Gut
3/8 circle Reverse Cutting
12/pkg \$31.95

Medium - 19mm needle

PS-2 / P12	3-0 18" long	875-2085
	3-0 27" long	875-2086
	4-0 18" long	875-2087

Medium - 16mm needle

PS-3 / P11	4-0 18" long	875-2088
	5-0 18" long	875-2089



SURGICAL

Sutures

nylon

(Applied Dental)

Unify Nylon 3/8 circle Reverse Cutting 18" long 12/pkg \$16.95

Medium - 19mm needle

PS-2/FS-2	3-0	875-2020
	4-0	875-2021
	5-0	875-2022



Medium - 16mm needle

PS-3	4-0	875-2023
	5-0	875-2024

(Applied Dental)

Unify Premium Nylon 3/8 circle Reverse Cutting 18" long 12/pkg \$29.95

Medium - 19mm needle

PS-2 Blk	3-0	875-2025
	4-0	875-2026
	5-0	875-2027



PS-2 Undyed

	3-0	875-2033
	4-0	875-2034
	5-0	875-2035

Medium - 16mm needle

PS-3 Blk	4-0	875-2036
	5-0	875-2037

PS-3 Undyed

	4-0	875-2038
	5-0	875-2039

Non-Absorbable, Nylon (Covidien)

12/pkg C-13 3-0 875-3100 \$28.95



Non-Absorbable, Nylon, Black Monofilament (Look)

3/8 Circle Reverse Cutting, 12/pkg

922B	C-6	4-0	875-4870	\$29.95
918B	C-6	5-0	875-4871	\$29.95
925B	C-7	3-0	875-4872	\$29.95
931B	C-17	4-0	875-4873	\$27.99
913B	C-17	5-0	875-4874	\$27.99
911B	C-17	6-0	875-4875	\$27.99
915B	C-22	6-0	875-4876	\$35.99



Sutures

chromic gut

(Sharp)

Absorbable, Chromic Gut \$19.95
3/8 Circle Reverse Cutting, 12/pkg

561B	C-7	3-0	875-5030
595B	C-6	4-0	875-5034
525B	C-31	4-0	875-5035



(Applied Dental)

Unify Premium Chromic Gut 3/8 circle Reverse Cutting 18" long 12/pkg \$29.95

Medium - 19mm needle

PS-2	3-0	875-2045
	4-0	875-2046



Medium - 16mm needle

PS-2	4-0	875-2047
	5-0	875-2048

Absorbable, Chromic Gut (Look)
3/8 Circle Reverse Cutting, 12/pkg \$37.95

560B	C-6	3-0	875-4880
559B	C-6	4-0	875-4881



SURGICAL

Sutures

chromic gut

Absorbable, Chromic Gut
(Covidien)



36/pkg	C-13	3-0	875-3000	\$69.95
36/pkg	C-23	3-0	875-3002	\$69.95

Sutures

PGA

(Applied Dental)

Absorbable, Unify PGA
Polyglycolic Acid, 3/8 circle Reverse
Cutting 18" long, 12/pkg



Medium - 19mm needle		\$31.95
PS-2/FS-2 Violet		
	3-0	875-2070
	4-0	875-2071
	5-0	875-2072
PS-2/FS-2 Undyed		
	3-0	875-2073
	4-0	875-2074
	5-0	875-2075

Sutures

polypropylene

(Applied Dental)

Unify Polypropylene
3/8 circle Reverse Cutting 18" long
12/pkg

\$25.95



Medium - 19mm needle		
PS-2/FS-2	3-0	875-2040
	4-0	875-2041
	5-0	875-2042

Medium - 16mm needle

PS-3	4-0	875-2043
	5-0	875-2044

(Applied Dental)

Unify Premium Polypropylene
3/8 circle Reverse Cutting 18" long
12/pkg

\$39.95



Medium - 19mm needle		
PS-2	4-0	875-2055
	5-0	875-2056

Medium - 16mm needle

PS-3	4-0	875-2057
	5-0	875-2058

Medium - 16mm needle

PS-3 Violet		
	4-0	875-2076
	5-0	875-2077
PS-3 Undyed		
	4-0	875-2078
	5-0	875-2079

(Applied Dental)

Absorbable, Unify Premium PGA
Polyglycolic Acid, 3/8 circle Reverse
Cutting 18" long, 12/pkg

\$39.95

Medium - 19mm needle

PS-2	3-0 Violet	875-2130
	4-0 Undyed	875-2131
	4-0 Violet	875-2132

Medium - 16mm needle

PS-3 Violet		
	4-0	875-2135
	5-0	875-2136
PS-3 Undyed		
	4-0	875-2137
	5-0	875-2138



Resorbable Tissue Replacement



Resorbs. Regenerates. Restores.



R.T.R.

Absorbable Bone Regeneration and Membrane



SURGICAL



Absorbable, Unify RPGA
Rapid Polyglycolic Acid 3/8 circle
Reverse Cutting 18" long, 12/pkg

Medium, 19mm		\$29.95
FS-2 / C-6	3-0 875-2090 4-0 875-2091	
Small, 13mm		
P-3 / C-3	4-0 875-2100 5-0 875-2101 6-0 875-2102	



Needle Chart

- C-3 13mm 3/8 circle
- C-6 18mm 3/8 circle
- C-7 24mm 3/8 circle
- C-8 26mm 3/8 circle
- C-17 12mm 3/8 circle
- C-22 16mm 3/8 circle
- C-26 15mm 1/2 circle
- C-31 24mm 1/2 circle



Sutures

PDO

(Applied Dental)

Absorbable, Unify PDO
Polydioxanone, 3/8 circle Reverse Cutting
18" long 12/pkg

Medium - 19mm needle		\$31.95
PS-2/FS-2 Violet	3-0 875-2110 4-0 875-2111 5-0 875-2112	



Medium - 16mm needle

PS-3 Violet	4-0 875-2115 5-0 875-2116
-------------	----------------------------------------

(Applied Dental)

Absorbable, Unify Premium PDO
Polydioxanone, 3/8 circle Reverse Cutting
18" long 12/pkg

Medium - 19mm needle		\$39.95
PS-2/FS-2 Violet	4-0 875-2125 5-0 875-2126	


Medium - 16mm needle

PS-3 Violet	4-0 875-2128 5-0 875-2129
-------------	----------------------------------------

Bone Graft Materials

allograft

InterOss (SigmaGraft)
100% Natural, Biocompatible,
Highly Purified, Large Inner
Surface Area, Long Term
Stability, Multi-Porosity



Vial	0.25-1 mm
0.25g	870-1125
0.5g	870-1126
1.0g	870-1127
2.0g	870-1128
Vial	1-2 mm
0.5g	870-1130
1.0g	870-1131
2.0g	870-1132

Syringe	0.25-1 mm	1-2 mm
0.125g	870-1135	-
0.25g	870-1136	870-1137



Bone Graft Material

allograft

Ti-oss Vial (Ti-oss)
made from 100% cancellous bone without any cortical portion. Innovative pulverizing technique allows multi-porous structure, maximizing blood vessel ingrowth



Vial	0.5-1.2 mm	1.2-1.7 mm	0.2-1.0 mm	
0.25g	870-1155	870-1160	870-1170	\$64.95
0.5g	870-1156	870-1161	870-1171	\$89.95
1.0g	870-1157	870-1162	870-1172	\$147.95
2.0g	870-1158	870-1163	870-1173	\$243.95

Ti-oss Syringe (Ti-oss)
made from 100% cancellous bone without any cortical portion. Innovative pulverizing technique allows multi-porous structure, maximizing blood vessel ingrowth



	0.5-1.2 mm	1.2-1.7 mm	
0.25g	870-1175	-	\$75.95
0.5g	870-1176	-	\$99.95
0.25g	-	870-1178	\$75.95
0.5g	-	870-1179	\$99.95

AlloOss Cortical Particulate (Ace Surgical)

Offers superior structural strength, which helps to maintain space and volume during the remodeling phase



Vial	0.25-1.0 mm	1.0-2.0 mm	
0.5cc	870-1190	870-1195	\$77.95
1.0cc	870-1191	870-1196	\$101.95
2.0cc	870-1192	870-1197	\$170.95

Bone Graft Materials

synthetic

R.T.R. Syringe (Septodont)
β-tricalcium phosphate bone grafting material



Curved Syringe
0.8 cc of β-tricalcium phosphate granules (diameter 0.5 to 1 mm) in sterile syringe (individually packaged).


870-1300 \$109.95

Unigraft® (Unicare)
Best alternative to Autogenous Bone, Resorbable. Unigraft is made of the highest quality bioactive glass by a proprietary process. The synthetic bioactive granules are optimally sized for ease of administration and enhancement of tissue growth. Bioactive Glass 200-420 Micron



Standard Kit
0.37g Vial x 5 pkg 870-1000 \$164.95
Bulk Kit
1g Vial x 5 pkg 870-1010 \$268.95

DentoGen (Orthogen)
Caclium sulfate hemyhydrate, post extraction, apioectomy & sinus lifts



Patient kits \$132.95
2 x 1 g of DentoGen material
Regular set 870-1040
Fast set 870-1045

NanoGen (Orthogen)
Nanocrystalline caclium sulfate hemyhydrate, post extraction, apioectomy & sinus lifts \$138.95



Patient kits
2 x 1 g of NanoGen material
Regular set 870-1050
Fast set 870-1055



SURGICAL

DM Bone (Meta Biomed)
Synthetic resorbable material
Biphasic calcium phosphate bioceramics
for implants



Vial	0.1-0.3 mm	0.3-0.5 mm	0.5-1.0 mm	1.0-2.0 mm	
0.25g	870-1250	870-1260	870-1270	870-1280	\$49.95
0.5g	870-1251	870-1261	870-1271	870-1281	\$59.95
1.0g	870-1252	870-1262	870-1272	870-1282	\$64.95
2.0g	870-1253	870-1263	870-1273	870-1283	\$69.95

Bone Graft Materials

membrane - Resorb

R.T.R. Membrane (Septodont)
Absorbable Collagen Membrane

R.T.R. Membrane is an absorbable, implantable material that is indicated for guided tissue regeneration procedure in periodontal defects to enhance regeneration of the periodontal apparatus.



15 x 20mm Absorbable Collagen Membrane	870-1301	\$79.95
20 x 30mm Absorbable Collagen Membrane	870-1302	\$93.95

(Unicare)

Cytoflex® Resorb
Resorbable Synthetic Barrier Membrane
Space-making barrier in guided tissue regeneration procedure. Provide a favorable environment for neo-vascularization and repopulation of bone cells to regenerate osseous tissue. The flexible barrier membranes are easily adaptable to tissue contours. Composed of synthetic polylactide and polyglycolide copolymers.



Single Pack		
20mm x 25mm	870-1150	\$83.95
30mm x 40mm	870-1151	\$109.95


HeliCote™ (Integra/Miltex)
Used in palatal donor sites or mucosal flaps
3/4" x 1 1/2" 870-1030 \$134.95
pack of 10



HeliTape™ (Integra/Miltex)
Used for dressing minor wounds. Can be used to aid in the closure of graft sites. Can also help repair the Schniderian Membrane. Dressing retains structure integrity even when wet. HeliTape® is highly porous and absorbs blood and wound exudate.



1" x 3"		
pack of 10	870-1020	\$187.95



R.T.R.
Absorbable Bone Regeneration and Membrane



Bone Graft Materials

membrane

HeliMend™ (Integra/Miltex)

Absorbed in 4 - 8 weeks

The HeliMend™ and HeliMend™ Advanced absorbable collagen membranes are made from purified bovine tendon. The product provides wound stabilization and creates space during guided tissue regeneration procedures where a longer absorption rate is desired.



Single Pack

1.5cm x 2cm	870-1200	\$86.95
2cm x 3cm	870-1201	\$106.95
3cm x 4cm	870-1202	\$173.95

HeliMend™ Advanced

Absorbed in 18 weeks

Single Pack

1.5cm x 2cm	870-1210	\$101.95
2cm x 3cm	870-1211	\$121.95
3cm x 3cm	870-1212	\$178.95



HeliPlug™ (Integra/Miltex)

Used to control bleeding and develop blood clots

3/4" x 3/4"	870-1035	\$121.95
-------------	----------	----------

pack of 10



GelFoam #4 (Upjohn)

Absorbable gelatin sponge 9039605 size 4, 2 x 2cm

6 envelopes x 2 sponges each	857-4500	\$254.95
------------------------------	----------	----------



Hemostatic Gauze (ActCel)

HG020022 Used for extractions and dry socket prevention, perio, endo and all oral surgical procedures 2" x 2"

20/box	857-4510	\$109.95
--------	----------	----------



Hemostatic Agents

BenaCel (Unicare)

Sterile, absorbent hemostatic gauze Stops bleeding quickly, prevents dry socket formation, completely resorbable, non-allergenic and non-immunogenic.

C-001 0.5cm x 0.7cm, Plug		
10/pack	857-1201	\$32.95
C-002 1.5cm x 1.5cm, Sheet		
12/pack	857-1202	\$44.95
C-003 0.6cm x 0.8cm, Plug		
8/pack	857-1203	\$39.95
C-005 5cm x 5cm, Sheet		
10/pack	857-1205	\$54.95



Splinting

Splintmat & Performat (Pulpdent)

1 meter/39", for permanent stabilization of teeth

SPL-C, Splintmat, coarse woven nylon mesh grid roll, 5.5mm	857-1455	\$64.95
------------------------------------------------------------	----------	---------

SPL-F, Splintmat, fine woven stainless steel mesh grid roll, 4mm	857-1460	\$64.95
------------------------------------------------------------------	----------	---------

SPL-P, Performat, perforated stainless steel strip roll, 3.5mm	857-1465	\$64.95
----------------------------------------------------------------	----------	---------



Fax your order to: 310.537.1502



SURGICAL

Irrigation/Injectable

bag/bottle

0.9% Sodium Chloride U.S.P.
(Baxter)

Plastic Bottle

1000 ml	857-1400	\$9.95
500 ml	857-1401	\$6.95



Saline/ Sterile water (Baxter)

Plastic Bottle

1000 ml	857-1415	\$9.95
---------	----------	--------



Saline/ Sterile water (Braun)

Plastic Bottle, Irrigation

1000 ml	857-1414	\$9.95
---------	----------	--------

0.9% Sodium Chloride U.S.P.
(Hospira)

Plastic Bottle

1000 ml	857-1405	\$11.95
500 ml	857-1406	\$9.95



5% Dextrose Solution

(Brown)

Plastic Bag

1000 ml	857-1425
500 ml	857-1424



0.9% Sodium Chloride Injectable

(Baxter)

Plastic Bag

1000 ml	857-1411	\$9.95
500 ml	857-1412	\$5.95



Medication

0.9% Sodium Chloride Injectable

(Hospira)

Plastic Bag

1000 ml	857-1410	\$9.95
500 ml	857-1409	\$9.95



Amoxicillin Capsule 500mg

(Caramel)

500/Pk	857-1450	\$69.95
--------	----------	---------



Cephalexin Capsules 500mg

(Ascend)

500/Pk	857-1451	\$79.95
--------	----------	---------



Saline/ Sterile water (Baxter)

Plastic Bag

1000 ml	857-1416	\$9.95
500 ml	857-1417	\$7.95



Doxycycline Hyclate 50/pk (Teva)

50ml	857-1460	\$49.95
100ml	857-1461	\$39.95



Saline/ Sterile water (Hospira)

Plastic Bag, Irrigation

1000 ml	857-1418	\$9.95
500 ml	857-1419	\$8.95

Ibuprofen Tablets 800mg

(Amneal) 500/Pk

400 mg	857-1456	\$39.95
600 mg	857-1455	\$49.95
800 mg	857-1453	\$59.95

