

Septalloy is ideal for:

- Any posterior restoration
- Core build-up for crown preparation
- Core build-up with pre-fabricated posts and pins

Septalloy — the zinc free amalgam restorative that's a viable and cost effective alternative to composites.

800-872-8305 | septodontusa.com

Septalloy is the single-composition, non-gamma II, dental amalgam with a track-record for high performance and durability. A ternary alloy with high copper single composition, Septalloy exhibits excellent compressive strength and maintains marginal integrity over time. Its reliable consistency and easy handling make Septalloy the perfect solution to the demands of large posterior restorations as well as core build-ups for crown preparation and pre-fabricated posts.



From the manufacturers of Septocaine®

Amalgams
capsules

(Carl Heyer)

Spherodon M Spherical

	1-Spill	2-Spill	3-Spill
50 Capsules	\$59.95	\$69.95	\$89.95
Reg	125-1630	125-1631	125-1632
Fast	125-1633	125-1634	125-1635
500 Capsules	\$399.95	\$579.95	\$637.95
Reg	125-1640	125-1641	125-1642
Fast	125-1643	125-1644	125-1645



Spherodon + Spherical

	1-Spill	2-Spill	3-Spill
50 Capsules	\$59.95	\$69.95	\$89.95
Reg	125-1650	125-1651	125-1652
Fast	125-1653	125-1654	125-1655
500 Capsules	\$499.95	\$699.95	\$799.95
Reg	125-1660	125-1661	125-1662
Fast	125-1663	125-1664	125-1665



Nogama Dispersed

	1-Spill	2-Spill	3-Spill
50 Capsules	\$59.95	\$69.95	\$92.95
Reg	125-1670	125-1671	125-1672
Fast	125-1673	125-1674	125-1675
500 Capsules	\$552.95	\$697.95	\$919.95
Reg	125-1680	125-1681	125-1682
Fast	125-1683	125-1684	125-1685



Dexiter (Carl Heyer)

Dispersed, Regular

	1-Spill	2-Spill	3-Spill
50 Capsules	\$59.95	\$79.95	\$99.95
	125-1606	125-1607	125-1608



Spherical, Regular

	1-Spill	2-Spill	3-Spill
50 Capsules	125-1616	125-1617	125-1618
500 Capsules	125-1621	125-1622	125-1623



KOBRE40 (Dharma)

Economical Non-Gamma 2 Admix Alloy

- High Copper Content
- Zinc Free
- High Strength and Positive Dimensional Change ensures longevity and patient satisfaction
- Optima Handling Characteristics

	1-Spill	2-Spill	3-Spill
50 Capsules	125-3000	125-3001	125-3002
500 Capsules	125-3003	125-3004	125-3005



WykAlloy (Wykle Research)

Economical alloy with high strength
Long lasting
Crunchy feeling while condensing
High luster restorations
Excellent marginal adaption
No Gamma2- No corrosion
Pre-dosed for a consistent mix.
Composition: 43.1% Ag, 30.8% Sn and 26.1% Cu

	1-Spill	2-Spill	3-Spill
50 Capsules	\$39.95	\$59.95	\$79.95
Reg 5min	125-9131	125-9132	125-9133
500 Capsules	\$399.95	\$599.95	\$699.95
Reg 5min	125-9141	125-9142	125-9143



Fax your order to: 310.537.1502

Prices are subject to change without notice.



Septalloy NG50 (Septodont)

Zinc-free amalgam, single composition, non-gamma 2, ternary dental amalgam, easy handling and greater corrosion resistance. Mixing time: 6 to 8 sec working time: 3.5 min, Regular set

	1-Spill	2-Spill
50 Capsules	\$69.95	\$85.95
	125-1000	125-1001
500 Capsules	\$589.95	\$749.95
	125-1002	125-1003



Dispersalloy (Cauk Dentsply)

Dispersed phase alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$165.95	\$205.95	\$239.95
Reg	125-1050	125-1051	125-1052
Fast	125-1053	125-1054	125-1055
500 Capsules	\$1139.95	\$1544.95	\$2009.95
Reg	125-1056	125-1057	125-1058
Fast	125-1059	125-1060	125-1061



GS-80 (SDI)

High Strength Amalgam, 40% Silver. A non-gamma 2 admix alloy. GS-80 has very high early tensile and compressive strength.

	1-Spill	2-Spill	3-Spill
50 Capsules	\$115.95	\$136.95	\$174.95
Reg	125-2751	125-2752	125-2753
Fast	125-2754	125-2755	125-2756
500 Capsules	\$779.95	\$959.95	\$1098.95
Reg	125-2801	125-2802	125-2803
Fast	125-2811	125-2822	125-2833



Unison (Cauk Dentsply)

Spherical alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$159.95	\$205.95	\$254.95
	125-1100	125-1101	125-1102
500 Capsules	\$1624.95	\$2065.95	\$2554.95
	125-1103	125-1104	125-1105



Lojic+ Capsules (SDI)

Dispersed alloy

Permite plunger colors:

1-Spill: Orange, 2-Spill: Cream, 3-Spill: Brown, 5-Spill: Bright pink

	1-Spill	2-Spill	3-Spill	5-Spill
50 Capsules	\$145.95	\$187.95	\$323.95	\$306.95
Regular	125-7700	125-7701	125-7702	125-7703
Fast	125-7704	125-7705	125-7706	125-7707
500 Capsules	\$1014.95	\$1190.95	\$1425.95	\$2669.95
Regular	125-7708	125-7709	125-7710	125-7711
Fast	125-7712	125-7713	125-7714	125-7715



Megalloy® EZ (Cauk Dentsply)

Spherical alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$189.95	\$218.95	\$264.95
	125-1125	125-1126	125-1127
250 Capsules	\$589.95	\$692.95	\$882.95
	125-1128	125-1129	125-1130



Permite™ (SDI)

Dispersed Phase Alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$132.95	\$171.95	\$219.95
Reg	125-7616	125-7617	125-7618
Fast	125-7619	125-7620	125-7621
500 Capsules	\$1017.95	\$1157.95	\$1411.95
Reg	125-7610	125-7611	125-7612
Fast	125-7613	125-7614	125-7615



Tytin (Kerr)

Spherical alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$199.95	\$239.95	\$279.95
Regular	125-8051	125-8052	125-8053
Fast	125-8054	125-8055	125-8056
500 Capsules	\$1439.95	\$1724.95	\$2067.95
Regular	125-8101	125-8102	125-8103
Fast	125-8111	125-8122	125-8133



Tytin FC (Kerr)

Spherical alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$208.95	\$249.95	\$299.95
Regular	125-8251	125-8252	125-8253
Fast	125-8254	125-8255	125-8256
500 Capsules	\$1354.95	\$1629.95	\$1955.95
Regular	125-8261	125-8262	125-8263
Fast	125-8264	125-8265	125-8266



Phasealloy (Wykle Research)

65.9% Silver, disposable phase alloy

- Multiple set times, low Hg vapor release
- Extended, smooth curve
- Superior marginal integrity
- No gamma2 - no corrosion
- Low creep
- High compressive strengths
- Hg safe capsules design



	1-Spill	2-Spill	3-Spill
100 Capsules	\$129.95	\$179.95	\$159.95
Reg 5min	125-9020	125-9021	125-9022
Fast 3min	125-9025	125-9026	125-9027
500 Capsules	\$649.95	\$699.95	\$1009.95
Reg 5min	125-9002	125-9003	125-9004
Fast 3min	125-9005	125-9006	125-9007

Contour (Kerr)

Dispersed alloy

	1-Spill	2-Spill	3-Spill
50 Capsules	\$177.95	\$213.95	\$256.95
Regular	125-0906	125-0907	125-0908
Fast	125-0916	125-0917	125-0918
500 Capsules	\$1149.95	\$1379.95	\$1654.95
Regular	125-0901	125-0902	125-0903
Fast	125-0911	125-0922	125-0933



Original D (Wykle Research)

69.5% Silver, depression system.

- Consistent mixes with controlled set times
- Firm resistance to condensing
- Low creep
- Superior marginal integrity

- No gamma2 - no corrosion
- High early/long term compressive strengths
- Hg safe capsule design



	1-Spill	2-Spill	3-Spill
500 Capsules	\$609.95	\$879.95	\$1119.95
Reg 5min	125-9101	125-9102	125-9103
Fast 3min	125-9104	125-9105	125-9106
100 Capsules	\$134.95	\$184.95	\$234.95
Reg 5min	125-9107	125-9108	125-9109
Fast 3min	125-9110	125-9111	125-9112
50 Capsules	\$74.95	\$99.95	\$124.95
Reg 5min	125-9115	125-9116	125-9117
Fast 3min	125-9118	125-9119	125-9120

Zenith Type-T Caps (DMG)

Spherical alloy, Fast Setting

	1-Spill	2-Spill	3-Spill
50 Capsules	\$67.95	\$82.95	\$104.95
	125-7750	125-7751	125-7752
500 Capsules	\$649.95	\$789.95	\$969.95
	125-7753	125-7754	125-7755



Zenith Premium Caps (DMG)

Dispersed phase alloy, Fast Setting

	1-Spill	2-Spill	3-Spill
50 Capsules	\$67.95	\$82.95	\$104.95
	125-7770	125-7771	125-7772
500 Capsules	\$649.95	\$789.95	\$964.95
	125-7773	125-7774	125-7775



Amalgamators

Amalgamator - Digital

First Choice DB338
 Quiet, Digital 110V Direct drive
 motor mixing system 135-3025 \$195.95



Ultramat2 (SDI)
 Quiet, high speed and efficient mix-
 ing system 135-7500 \$709.95



Optimix (kerr)
 3,000 - 4,800 Cpm speed range, quiet oper-
 ation, 3-touch system, quiet operation
 135-2000 \$1349.95



(Wykle Research)
 Consistent Mixes with low mix times, superior
 amalgam properties, indirect drive motor,
 external fuse block, four-year warranty
 Max Ultra High Speed Mixer
 135-3000 \$429.95



A-1000 Duo-Amalgamator (Vector)
 Dual speed digital amalgamator
 auto-shut off door, quiet
 135-3015 \$295.95



Amalgam Carrier

(FCI)
 Double End \$7.95
 Medium \varnothing 0.2mm / Small, 1.5 mm 130-1031
 Medium \varnothing 0.2mm / Large, 2.5 mm 130-1032
 Medium \varnothing 0.2mm / Jumbo, 3.0 mm 130-1033



(Arnel)
 Double End, Stainless Steel \$14.95
 Reg/large 130-0050



(Miltex)
 Double-End carrier \$46.95
 Stainless Steel Nylon
 Med-Large 130-5100 130-5104
 Med-Jumbo 130-5101 130-5105
 Large-Jumbo 130-5102 130-5106
 Jumbo-Jumbo 130-5103 130-5107



Single-End carrier \$28.95
 Stainless Steel Nylon
 Medium 130-5108 130-5110
 Jumbo 130-5109 130-5111



SCT Carbide Tip Amalgam Carrier (Miltex)


Double-End carrier \$83.95
 Med-Large 130-5150
 Med-Jumbo 130-5151
 Large-Jumbo 130-5152
 Jumbo-Jumbo 130-5153
 Single-End carrier \$49.95
 Medium 130-5154
 Large 130-5155
 Jumbo 130-5156



(Nordent) \$53.95
 Unconditional 2 year guarantee,
 double end, stainless steel
 Regular/Mini 130-5010
 Regular/Large 130-5011
 Regular/Jumbo 130-5012
 Jumbo/Jumbo 130-5013



Neverclog™ Amalgam Gun (Premier)
 Plastic, autoclavable, self lubricating \$28.95
 Mini right angle, 1.25mm 130-5180
 Right angle, 2.3mm 130-5181
 Straight 2.3mm 130-5182



Amalgam Gun (Accent on Dental)
 Plastic, autoclavable, self lubricating \$14.95
 Mini right angle, 1.25mm 130-5190
 Right angle, 2.3mm 130-5191
 Straight 2.3mm 130-5192



Amalgam Squeeze Cloth

(JR Rand)
 Made of fine, closely woven cotton
 Round 3" diameter
 500/pkg 140-0000 \$10.95



Amalgam Wells

(FCI)
 Stainless Steel, non-slip silicon "O" ring base
 145-1030 \$5.95



(Arnel)
 Standard #1 145-0000 \$6.50
 Stainless Steel, has a non-slip silicon "O"
 ring base
 Non-Skid #2 145-0100 \$7.50
 With Diamond Knurl & Rubber Base
 Mini Non-Skid #3
 With Rubber Base 145-0200 \$7.50



(Miltex)
 Non-slip 145-0250 \$19.95
 Small
 fits Thompson Cassettes 145-0251

