


CLEARFIL™ SA CEMENT – a new self-adhesive composite resin cement

Structure of the adhesive monomer MDP



The well-known and proven monomer MDP, 10-Methacryloyloxydecyl dihydrogen phosphate, is used for the bond strength, which displayed very good adhesion to hard tooth substance as well as metal and zirconium oxide in a number of trials.

Kuraray Medical Inc. was able to prove this high tolerance towards moisture for the bonding on enamel and dentin.



CEMENT & LINERS

Temporary Cement

powder & liquid

Hy-Bond™ (Shofu)

Polycarboxylate temporary cement

Standard Kit PN 1176 240-0400 **\$79.95**

60g powder, 34ml liquid, dispensers & instructions

Bulk Kit PN 1173 240-0405 **\$109.95**

125g powder, 60ml liquid, dispensers & instructions



Hy-Bond™ Zoe (Shofu)

Zinc Oxide eugenol temporary filling

Kit #PN 1180 240-0410 **\$29.95**

20g powder, 12ml liquid, dispenser & instructions

Powder #PN 1181, 20g 240-0415 **\$19.95**

Liquid # PN 1182, 12ml 240-0416 **\$12.95**



tube

Temp-Bond (Kerr) **\$34.95**

50g base, 15g accelerator, mixing pad

Metal tube 240-0500

Plastic tube 240-0501



Temp-Bond NE (Kerr) **\$34.95**

Non-eugenol

50g base, 15g accelerator, mixing pad

Metal tube 240-0510

Plastic tube 240-0511



Zone Temporary Cement (DUX)

Standard Package 27040

Base, Catalyst and mixing pad

240-0600 **\$36.95**



RelyX™ Temp NE (3M ESPE)

Temporary cement Zinc Oxide

Standard Pkg

36g base & catalyst plate, mixing pad

Eugenol 56660 240-3620 **\$28.95**

Non-Eugenol 35015 240-3621 **\$35.95**



unit dose

Embonte (DUX)

ZOE / long term cementation of temporary crowns splints & trials, 30 sec mix time and working time of 2 min

Unit Dose #27202 240-1021 **\$19.95**

25 x 1.5g pouches



Zone - Unit Dose (DUX)

Temporary Cement NE 27039

25 x 1.5g pouches 240-0800 **\$19.95**



Temp-Bond Unidose (Kerr)

Perfect base/accelerator mix ratio each time, which minimizes material waste and maximizes convenience

50/pkg 240-1071 **\$34.95**



Temp-Bond NE Unidose (Kerr)

non-eugenol

50/pkg 240-1072 **\$34.95**



automix/syringe

PacTemp Automix NE (Pac-Dent)

Non-eugenol resin based temporary cement

Kit 240-0700 **\$19.50**

1 x 6ml dual barrel syringe & 10 mixing tips



PacTemp NE (Pac-Dent)

Non-eugenol resin based temporary cement, comes in separate syringes

Kit 240-0710 **\$32.50**

10ml Part A & 10ml Part B syringe



Embonte Syringe (DUX)

ZOE / long term cementation of temporary crowns splints & trials, 30 sec mix time and working time of 2 min

#27200 15g syringe 240-1020 **\$36.95**



Embonte Automix Syringe (DUX)

ZOE / long term cementation of temporary crowns splints & trials, 30 sec mix time and working time of 2 min

#27203 4 grams, 8 tips 240-1025 **\$36.95**



Zone Automix Syringe (DUX)

Temporary Cement Automix
Non-eugenol, 4g syringe & 8 tips **\$36.95**

Shade A1 27046 240-1040
Universal 27047 240-1041



Zone Temporary Cement (DUX) Syringe Delivery

30 sec mixing time & 2 min working
time, Non-eugenol, 15g **\$39.95**

Universal 27043 240-1000
Shade A1 27045 240-1001



TNE Dual Syringe (Temrex)

Non-eugenol temporary cement

Standard Pack 240-1050 **\$32.95**

6g dual-barrel syringe each parts
A & B, 15 tips, & releasing agent



Nextemp (Premier)

Non-eugenol resin-matrix formulation to
combine the multiple benefits of fluoride
release to protect tooth structure

Standard Kit **\$24.95**

5ml Automix dual barrel syringe &
10 mix tips 240-0900

Bulk Pack **\$49.95**

4 x 5ml Automix dual
barrel syringe 240-0910



Integrity® Tempgrip (Dentsply)

Temporary Crown & Bridge cement. Zinc
Oxide, non-eugenol

Standard Pkg 240-3550 **\$79.95**
2 x 9g syringes & 20 mixing tips

Mixing Tips, 50/pkg 240-3551 **\$54.95**



PowerTemp (J Morita)

Standard Pack 240-1090 **\$59.95**
50 ml cartridge & 10 blue small mix tips



Sensitemp Resin (Sultan)

Regular resin temporary cement
4 ml, w/20 tips 240-1060 **\$29.95**

Sensitemp Nezo (Sultan)

Non-eugenol zinc oxide temporary cement
4 ml w/tip 240-1061 **\$29.95**

Sensitemp Zoe (Sultan)

Eugenol zinc oxide temporary cement
4 ml w/tip 240-1062 **\$29.95**



Temp-Bond Original (Kerr)

Automix, 2 syringes (11.8g each)
20 Mixing Tips 240-1070 **\$64.95**



Temp-Bond NE (Kerr)

Non-eugenol, automix
2 x 11.8g syringes & 20 Mixing Tips
240-1075 **\$64.95**



Temp-Bond Clear (Kerr)

Eugenol-free, dual-cure, fluoride-
releasing temporary resin, automix
7g syringe, 10 Mixing Tips
240-1080 **\$49.95**



TempoSil (Coltene Whaledent)

Unique addition cured, silicone based, zinc
oxide temporary cement with silane for
optimal adhesion and superior marginal
integrity. Direct, easy application
Time saving, unique silicone formula, can
be hand mixed for 20sec

Intro Kit **\$37.95**
2 x 5ml & tips 240-3575

Standard Pkg **\$59.95**
4 x 5ml 240-3576

Mixing Tips **\$29.95**
40/bag w/oral tip 240-3577



CEMENT & LINERS

Temporary Cement

automix/syringe

TempoSil®2 (Coltene Whaledent)
Temporary cement w silane agent for optimal adhesion & superior marginal integrity, automix

Trial Kit \$34.95

5ml each of Dentin & White, & 15 mixing tips w attached oral tips 240-3600

Intro Kit \$34.95

2 x 5ml each of Dentin & White, & 15 mixing tips w attached oral tips

Dentin 240-3601

White 240-3602

Refills, 4 x 5ml syringes \$59.95

Dentin 240-3603

White 240-3604



TurboTemp 2 Cement (Danville)

Non-eugenol, provides optimal cementation for provisional & permanent crowns

Kit 240-2085 \$26.95

10 ml automix syringe, & 6 waste saver mix tips



Temporary Resin Cement (Mizzy)

Low irritation potential, easily removed, optimize morshy time, 0.2°

2 Automix syringes 240-2020 \$39.95



Zoe Temporary Cement (Mizzy)

2 Automix syringes

Neutral color, no effect on aesthetics

Pulp protection from shock, cavity liner & temporary filling 240-2000 **\$49.95**



Zinc Oxide Temporary Cement (Mizzy)

Non-eugenol formula

Excellent marginal seals & adhesion

2 Automix syringes 240-2010 \$49.95



Temporary Filling

regular

Cavit Original (3M ESPE)

Light Pink

4 x 7g tubes 240-2040 **\$18.95**

28g Jar 240-2041 **\$18.95**



Cavit G (3M ESPE)

Gray and soften setting form

28g jar 240-2030 **\$18.95**



Cavit W (3M ESPE)

White and medium setting form

28g jar 240-2042 **\$18.95**



MD-Temp (Meta)

40g jar

\$4.95

Pink 240-2050

White 240-2051



IRM (Caulk)

Zoe intermediate restorative

Standard Pack 240-2060 \$49.50

38g powder + 14ml liquid

14ml liquid 240-2061 **\$19.50**

38g powder 240-2062 **\$34.95**



Zoe B&T (Caulk)

Zinc oxide eugenol base

25g powder, 14ml liquid

240-2080 **\$59.95**



unit dose

Tempit® Unidose (Centrix)
Single dose syringe tips for easy placement, maximum asepsis

Standard pkg 240-2021 **\$39.95**
40 x 0.35g tips



Tempit® Ultra-F (Centrix) **\$69.95**

Microfilled pylover resin with fluoride release. 30 x 0.25g pre-filled tips

Regular pkg 240-2022
w/Fluoride 240-2023



IRM Caps (Dentsply Caulk)
To activate/mix use amalgamator

50 x 0.34g capsules 240-2070 **\$79.95**



Temporary Filling

light-cure

Spacer™ (Vericom) **\$29.50**

Clear (Yellow) 240-2090
Blue 240-2091



Fermit™ (Vivadent) **\$49.50**

Temporary filling after inlay preparations
3 x 2.5g syringes 240-3000



Fermit-N™ (Vivadent) **\$49.50**

Non-eugenol
3 x 2.5g syringes 240-3001

First Fill® (Pentron)

Light cured temporary filling material, remains elastic after curing to facilitate removal without damaging the margins

Kit 240-3020 **\$39.95**
2 x 2ml syringes, instructions & MSDS



Tempit® L/C (Centrix)
Natural, translucent color

30 x 0.25g pre-filled tips 240-3660 **\$69.95**



Soleil™ (DMG)

Temporary resin material universal shade
2 x 2.5 g syringes 120-2100 **\$44.95**



Dental Stopping

Temporary Stopping (Diadent)

Use as a temporary filling

	30g	120g
White	246-4100	246-4104
Yellow	246-4101	246-4105
Pink	246-4102	246-4106
Assorted	246-4103	246-4107



Hygienic Temporary Filling

(Coltene/Whaledent)

4oz/pack **\$19.95**

White	246-4110
Pink	246-4112
Beige	246-4114



Permanent Cement

powder & liquid

Siment-Z (Deepak)

Zinc phosphate cement

Standard Kit 246-4050 **\$5.49**

32g powder, 17.5ml liquid

Powder, 32g 246-4051 **\$3.95**

Liquid, 17.5ml 246-4052 **\$2.95**



Durelon Bottle Kit (3M ESPE)

Carboxylate cement, chemically adheres to enamel, dentin, stainless steel

Triple Size Pack 242-1610 **\$129.95**

60g powder, 40g liquid, mix pad

Triple Liquid, 40g 242-1620 **\$79.95**

Triple Powder, 60g 242-1630 **\$59.95**



Permanent Cement

powder & liquid

RelyX Luting Cement (3M ESPE)

Glass Ionomer luting cement
for crown & bridge

Kit 242-1700 **\$129.95**

16g powder, 9ml liquid, & mixing pad



M-Bond (Tokuyama)

Self-etching adhesive resin cement

Kit 242-2000 **\$199.50**

Clear Powder 3g, Ivory 3g, Liquid 8ml,
Primer A & B 3ml each & accessories

Refill, 3g powder **\$19.95**

Clear 242-2010

Ivory 242-2011

Liquid, 8cc 242-2012 **\$79.95**



Fleck's (Mizzy)

Zinc oxyphosphate cement, for cementing
crowns, inlays, fixed bridge-work, orthodontic
fans, temporary fillings and as a base
under any filling material

Liquid

15cc 242-3000 **\$9.95**

35cc 242-3001 **\$17.95**

Powder, 29g **\$8.95**

Golden Yellow 242-3002

Light Yellow 242-3003

Yellow 242-3004

Snow White 242-3005



Ketac-Cem Radiopaque (3M ESPE)

Glass ionomer fluoride releasing luting
cement

Intro Pack 242-3100 **\$99.50**

33g powder + 12ml liquid

Triple Pack 242-3101 **\$199.50**

99g powder & 36ml liquid

12ml liquid 242-3102 **\$19.95**

33g powder 242-3103 **\$79.95**



Zinc Phosphate Cement (Dentonics)

Standard Pack **\$8.95**

32g powder & 17.5ml liquid

Light Yellow 242-0950

White 242-0951



GC Fuji™ I Final Cement C&B

(GC America)


For crown & bridge, enhanced self-cured
glass ionomer luting cement

1-1 pack 242-3150 **\$99.95**

35g powder + 25g liquid

Powder, 35g 242-3151 **\$69.95**

Liquid, 25g 242-3152 **\$29.95**



Fynal (Dentsply Caulk)

Zoe Cement

Standard Pack 242-3350 **\$69.95**

32g powder and 15 ml liquid

Powder 32g 242-3351 **\$49.95**

Liquid 15 ml 242-3352 **\$29.95**



C&B Metabond Quick (Parkell)

Complete Kit

Contains enamel etch, dentin activator,
base, catalyst, clear powder, tooth-
shade powder, ceramic mix shade,
brushes 242-3200 **\$295.50**



Hy-Bond™ (Shofu)

Polycarboxylate luting cement

Standard Kit PN 1160 242-3280 **\$79.95**

60g powder, 34ml liquid, dispensers
& instructions

Powder PN 1161, 60g 242-3281 **\$39.95**

Liquid PN 1162, 34ml 242-3282 **\$39.95**



Hy-Bond™ (Shofu)

Zinc Phosphate cement

Standard Kit PN 1170 242-3285 **\$29.95**

60g powder, 25ml liquid, dispensers
& instructions

Powder PN 1171, 60g 242-3286 **\$19.95**

Liquid PN 1172, 25ml 242-3287 **\$12.95**



Siment-P (Deepak)

Polycarboxylate cement

Standard Kit 242-3290 **\$6.95**

25g powder, 15ml liquid

Powder, 25g 242-3291 **\$3.95**

Liquid, 15ml 242-3292 **\$2.95**



GC Fuji Plus (GC America)

Resin reinforced multipurpose glass ionomer

1:1 Pack #431011 242-3161 **\$139.95**
15g powder, 7ml liquid, 6.5ml conditioner,
mixing pad & scoop

No-Etch Pack #431101 242-3160 **\$119.95**
as above but without conditioner

Powder, 15g 242-3170 **\$64.95**

Liquid, 7ml 242-3171 **\$49.95**

Conditioner, 6.5ml 242-3172 **\$37.95**



Permanent Cement

self-cure

Durelon Syringe Kit (3M ESPE)

Carboxylate cement

60g Powder, 12 syringe liquid CD,
mixing pad 242-3400 **\$99.95**

Triple Syringe 242-3410 **\$69.50**
12 x 40ml

Triple Powder 242-1630 **\$59.95**
60g jar



Embrace Wetbond Cement (Pulpdent)

Resin cement, automix syringe

Low Viscosity Kit **\$79.95**

7g syringe, 20 tips, 1.2ml each Seal-n-Shine
and OxOut 242-3450

Refill, 7g syringe 242-3451 **\$59.50**

High Viscosity Kit **\$79.95**

7g syringe, 20 tips, 1.2ml each Seal-n-Shine
and OxOut 242-3452

Refill, 7g syringe 242-3453 **\$59.50**

Accessories Kit 242-3454 **\$29.95**

Standard syringe, 1.2ml each of Seal-n-
Shine & OxOut with tips



Fuji Plus Capsule (GC America)

No-etch, resin reinforced multipurpose
glass ionomer

Starter Package 242-3505 **\$199.95**

50 capsules & capsule Applier III

Refill pack 242-3510 **\$169.95**

50 capsules, yellow shade

Capsule Applier III #135891

242-3515 **\$69.95**



Fuji I Capsules (GC America)

Light Yellow

Package of 50 242-3520 **\$169.95**



FujiCEM™ Automix (GC America)

The first paste-paste automix delivery
system, resin reinforced

Starter Kit #431573 240-3700 **\$219.95**

13.3g x 2-paste pack, 44 mixing tips,
& dispenser

Refill #2589 240-3701 **\$199.95**

13.3g x 2-paste pack, with mixing tips

Refill #2636 240-3702 **\$169.95**

13.3g x 2-paste pack, no mixing tips

Paste Pack Dispenser 240-3705 **\$59.95**



RelyX™ Unicem Maxicap (3M ESPE)

Ideally suited for multiple unit bridges.

Hygienic unidose (936g/capsule) consist-
ent mix w triturator, longer working
time - 2.30 min

Intro Kit 242-3530 **\$199.50**

20 Maxicap capsules, 10x A2 Universal, 6
translucent, 4 x A3 Opaque, activator &
applier, instrument set

Refills pkg of 20

Assorted 242-3540 **\$169.50**

10 x A2 Universal, 6 translucent, 4 x A3 Opaque

Translucent (20) 242-3550

A3 Opaque (20) 242-3551

Universal (20) 242-3552

Dispenser Maxicap System

Capsule applier 242-3555 **\$129.50**



RelyX™ Arc (3M ESPE)

Adhesive resin cement, for indirect
porcelain, ceramic, & cosmetic
restorations

Intro Kit 242-3560 **\$479.50**

Refill **\$99.50**

4.5g paste shade

A1 (Transparent) 242-3561

A3 (Universal) 242-3562

Cleaning Pads, 60/pack 242-3264 **\$9.95**



RelyX™ Unicem 2 Clicker Dispenser

(3M ESPE)

Self-adhesive universal resin cement in con-
venient clicker dispenser

Kit **\$199.50**

Transparent 242-3670

A2 Universal 242-3671

A3 Opaque 242-3672



Permanent Cement

self-cure

RelyX™ Luting Plus (3M ESPE)

Resin-modified Glass Ionomer luting cement for crown & bridge

Trial Kit 242-3500 **\$99.95**
1 x 11g syringe + Clicker dispenser

Refill Kit 242-3501 **\$199.95**
2 x 11g syringe



RelyX™ Unicem 2 Automix (3M ESPE)

Standard Refills **\$117.95**
8.5g (5ml) automix syringe, 15 mixing tips, 5 wide mixing tips, & 5 endo tips

Translucent 242-3700
A2 Universal 242-3701
A3 Opaque 242-3702



RelyX™ Unicem Aplicap (3M ESPE)

Universal cement product for metal & all-ceramic composite restorations

Intro Kit 56814 242-3580 **\$299.95**
50 aplicap capsules, 24 x A2 Universal, 16 translucent, 10 x A3 Opaque

Trial Kit 56839 242-3581 **\$89.50**
20 aplicap capsules, 10 x A2 Universal, 6 translucent, 4 x A3 Opaque, Applicap activator & applier set



Refills

50 Capsules: 295ml capsules **\$199.95**
Translucent 242-3585
A3 Opaque 242-3586
Universal 242-3587
Assorted 242-3588

Activator 242-3589 **\$59.95**

Applier 242-3590 **\$79.95**



Maxcem Elite™ (Kerr)

Self-Etch, self-adhesive cement

Intro Kit 242-4130 **\$139.95**
3 dual syringes: clear, white, yellow, 20 regular & 10 wide automix tips, 10 intra-oral tips, 10 root canal tips, instructions

Refills **\$99.95**

Syringes - 2 x 5g dual syringes, 20 tips, instructions

Clear 242-4131
White 242-4132
Brown 242-4133
White Opaque 242-4134
Yellow 242-4135



Smartcem2™ (Dentsply)

Self-adhesive cement w combination of strength & convenience

Syringe Intro Kit 242-3675 **\$197.95**

One of each shade: Light, Medium, Dark, Translucent, Opaque 5g, & 50 mixing tips

Refills **\$97.50**

Light 242-3676
Medium 242-3677
Dark 242-3678
Translucent 242-3679
Opaque 242-3680



PacCem Plus™ (Pac-Dent)

Designed primarily for use as a cavity liner and cement for orthodontic bands and dental prostheses

Kit 242-4260 **\$49.95**

6g syringe, 15 blue mixing tips

30 short blue mixing tips, 1:1

242-4261 **\$14.95**



Monocem™ (Shofu)

High retention values, eliminates sensitivity, high fluoride release, moisture tolerant and easy to clean-up, dual cure, with ultimate working time 100% polymerization and comes in an automix syringe. 12 micron thickness cement

Single Kit **\$79.50**

7ml syringe, 15 tips, 2 polishers

Bleach White 242-4150

Translucent 242-4151



SpeedCEM (Vivadent)

Self-adhesive self-curing universal resin cement w light-curing option

Kit **\$129.50**

9g syringe w 15 tips

Transparent 242-3880

Yellow 242-3881

Opaque 242-3882



Adherence® (Septodont)

Permanent resin cement with Fluoride

Kit 242-4175 **\$129.95**

2 x 4.5g automix syringes (Base & Catalyst), 10 mixing tips, mixing pad



Permanent Cement

dual-cure

Adherence® Dual Cure (Septodont)

Permanent resin cement with Fluoride

Kit 242-5390 **\$129.95**

2 x 4.5g automix syringes in Clear, White, Yellow, Brown, & Opaque shades, 25 mixing tips

Refills, 4g automix syringes \$34.95

- Clear 242-5391
- Yellow 242-5392
- Brown 242-5393
- White 242-5394
- Opaque 242-5395



Aegis C&B Cement (Bosworth)

Two-paste, dual-curing modified resin polymer material, with Amorphous calcium phosphate filler.

Standard Kit 242-5100 **\$169.95**

2 x 3ml syringe each part A & part B



Bistite® II DC (Tokuyama)

Self-etching dual-cure adhesive resin cement

Kit **\$159.50**

A & B paste 3.5g each, 5ml A & B primer 1, 5ml primer 2, air barrier, metaltite & accessories

- Clear 242-3600
- Ivory 242-3601
- Dentin 242-3602



Refill **\$69.50**

A & B paste 3.5g each

- Clear 242-3603
- Ivory 242-3604
- Dentin 242-3605

Refill, Metaltite

5ml bottle 242-3609 **\$49.95**

Breeze™ (Pentron)

Automix, single-step, dual-cure/light-cure or self-cure

Intro Kit 242-3790 **\$249.95**

3 x 4 ml syringe A2, 4 ml syringe Translucent 4 ml Opaceous White, 3 ml bottle of Silane, 40 automix tips, 10 intra-oral/automix tips, instructions & MSDS

Refills **\$59.95**

- A2 242-3791
- Translucent 242-3792
- Opaceous White 242-3793



Clearfil™ Esthetic Cement EX (Kuraray)

High-performance, esthetic, dual-cure composite cement system that includes a self-etch, single-step, dual-cure adhesive



Standard Kit \$199.50

Available in 3 shades, Double syringe for 1 shade, paste A & B, try-in paste for 1 shade, DC bond liquid A & B, ceramic primer, K-etchant gel, accessories

- Universal 2901KA 242-3710
- Clear 2903KA 242-3711
- Opaque White 2904KA 242-3712

Mini Kit \$149.95

1 x Double-syringe for 1 shade (Universal, Clear, and Opaque White): Paste A 2.4g, Paste B 4.4g, 1 x Try-in Paste for 1 shade 2.4g, 1 x DC Bond (Liquid A) 1ml, 1 x DC Bond (Liquid B) 1ml, 1 x Ceramic Primer 2ml, & accessories

- Universal #2911KA 242-3715
- Clear #2913KA 242-3716
- Opaque White #2914KA 242-3717

Refill - Paste \$99.95

each refill shade comes with Paste A 4.6g (2.3ml), Paste B 4.6g (2.3ml), & 20 mixing tips

- Universal #2921KA 242-3720
- Chroma #2922KA 242-3721
- Clear #2923KA 242-3722
- Opaque #2924KA 242-3723

Try-in paste \$29.95

Syringe 2.4g (1.8ml)

- Universal 2931KA 242-3725
- Chroma 2932KA 242-3726
- Clear 2933KA 242-3727
- Opaque 2934KA 242-3728

DC Bond Liquid A - 4ml 275-2521 **\$79.55**

DC Bond Liquid B - 4ml 275-2522 **\$79.55**

Ceramic primer, 4ml 242-3705 **\$69.95**

K-etchant gel, 6ml 288-2211 **\$9.95**

Clearfil Sa Cement™ (Kuraray)

self-adhesive resin cement, auto-mix delivery & fluoride-releasing. Does not require phosphoric acid treatment or any special primer or adhesive



Standard Kit - 8.3g syringe & 20 mixing tips \$79.50

- #2800KA Universal A2 242-3730
- #2801KA White 242-3731

Value Kit 3 x 4.6 ml (8.3g) syringes, 40 mixing tips, and 10 Endo tips \$199.50

- #2810KA Universal A2 242-3732
- #2811KA White 242-3733

Refill Tips \$12.50

- #2820KA 20 mixing tips 242-3734
- #2823KA 20 endo tips 242-3735

Permanent Cement

dual-cure

Cement-It™ w/Bond-1® C&B System

(Pentron)

Universal C&B Resin cement, Dual-cure

Kit 242-3794 **\$199.95**

4ml syringe A2, 4 ml Opaque White, 6ml Bond-1 primer/adhesive, 5ml PrepEze Desensitizer, 3ml syringe 37% etch gel, 3ml Silane, instructions & MSDS

Refills **\$59.95**

Automix syringe catalyst & base, 7g

A2 242-3801

Traslucent 242-3802

Opaque White 242-3803



G-Cem™ Capsule (GC America)

Self-adhesive resin luting cement in capsules, dual-cure

Starter Kit - Capsule 242-5350 **\$229.95**

50/box assorted capsule with dispenser

Refills **\$199.95**

50 capsules

A2 242-5351

AO3 242-5352

BO1 242-5353

Traslucent 242-5354

Assorted - Includes 20 x A2, 20 x Traslucent, 5 x AO3, 5 x BO1, instructions 242-5355

Capsules Applier 242-3515 **\$69.95**



G-Cem™ (GC America)

Self-adhesive resin luting cement dual-cure

Automix Kit **\$99.95**

2 x 4.8g Syringe (2.7ml), 20 x G-CEM Automix tips

Traslucent 242-5360

A2 242-5365

Regular Automix Tip

20 automix tips 242-3872 **\$24.95**



Geristore® Syringeable (DenMat)

The most biocompatible material for soft tissue & teeth. Dual-cure hydrophilic bis-GMA

Intro Kit 242-5400 **\$69.95**

5g syringe, 6 automixing tips, 6 intraoral tips, 3 curing rings, rule of 55 charts, & instruction cards



Infinity® Syringeable (DenMat)

Fluoride releasing, resin ionomer biocompatible dual cure cement. Its high bonding strength, tough fracture resistance, low film thickness, and radiopacity make Infinity the ideal materials for cementing PFM crowns and all gold restorations



Intro Kit #039690110 242-3746 **\$69.95**

5g syringe, 3 Auto-Mix Dispensing Tips, 3 curing rings, & accessories

Refill Kit #039670100 242-3748 **\$129.95**

10g syringe, 20 Auto-Mix Dispensing Tips 2 Intra-Oral Tips

Value Kit #039680100 242-3749 **\$329.95**

3 x 10g syringe, 60 Auto-Mix Dispensing Tips 4 Intra-Oral Tips

Panavia F2.0 (Kuraray)

Fluoride releasing and dual-cured for all ceramic restorations. The only fluoride releasing, anti-bacterial, self-adhering, self etching product on the market.

Complete Kit **\$249.95**

Available in 4 shades, 2.3ml each of Paste A & Paste B, 4ml, Primer II Liquid A & B, 6ml Oxiguard™ II, & accessories

TC-Tooth Color 501KA 242-5055

OP-Opaque 503KA 242-5056

EX-White 502KA 242-5057

Light 504KA 242-5058

Intro Kit **\$114.95**

Available in 4 shades, 1ml each of Paste A, Paste B, Primer II, Liquid A & B, 6ml Oxiguard II, & accessories

TC-Tooth Color 480KA 242-5050

OP-Opaque 482KA 242-5051

EX-White 481KA 242-5052

Light 483KA 242-5053

Refills

Paste A, 5g/2.3ml 242-5060 **\$49.95**

Paste B, 4.6g/2.3ml: **\$49.95**

Paste B

TC 242-5067

White 242-5068

Opaque 242-5069

Light 242-5070

ED Primer II, 4ml: **\$29.95**

Liquid A 242-5064

Liquid B 242-5065

6ml Oxiguard™ 242-5066 **\$24.95**



Prices are subject to change without notice.

Panavia 21 (Kuraray)

Allows one-step direct bonding to dentin and cut enamel without acid etching. Consistent mix with exact working and setting times.

Standard Kit **\$249.50**

Self-Cure Adhesive

7.9g catalyst, 7.6g Universal, 4ml ED Primer Liquid A & B, 6ml Oxyguard, 5ml etching gel, accessories

TC-Tooth Color	242-5000
OP-Opaque	242-5001
EX-White	242-5002

Refills **\$129.95**

7.9g Catalyst & 7.6g Universal paste

TC-Tooth Color	242-5003
OP-Opaque	242-5004
EX-White	242-5006

ED Primer **\$27.95**

Liquid A	242-5007
Liquid B	242-5008

Oxyguard™ II 6ml 242-5009 **\$24.95**

5ml etching agent 242-5011 **\$16.95**



Resilute (Puldent)

Kit

\$49.50

4g Base & 4g Catalyst

Universal Shade 242-5200

Clear Shade 242-5201



Variolink II (Ivoclar Vivadent)

Universal Dual Curing Cement

Standard Pack 242-5010 **\$199.50**

Contains 5 syringes of Base, 4 syringes of catalyst, and 1 mixing pad



Refills - 1 x 2.5g

Base **\$39.95**

Catalyst **\$39.95**

Bleach/XL-010 242-5020

Transparent 242-5021

White/100 242-5022

Opaque/100/White 242-5023

Universal/200/Yellow 242-5024

High Viscosity 242-5030

Low Viscosity 242-5031

High Viscosity Transp. 242-5032

Low Viscosity Transp. 242-5033

Lute-It® (Pentron)

Dual cure, used for cementing all ceramic & all composite indirect restorations

Kit 242-3890 **\$399.95**

8 x 3ml syringes Base shades & Try-In paste
2 x 3ml Catalyst, 1 of each 1ml syringe,
Etch gel, Bond-1, Silane, accessories

Refills - 3ml syringe

Base **\$29.95**

Try-In Paste **\$14.95**

	Base	Try-In Paste
Clear	242-3891	242-3900
Light	242-3892	242-3901
Dark	242-3893	242-3902
Opaque White	242-3894	242-3903
Universal	242-3895	242-3904
Lt Yellow	242-3896	242-3905
Dk Yellow	242-3897	242-3906
Gray	242-3898	242-3907
Opaque Universal	242-3899	242-3908



Implant Cement

Multilink Implant (Ivoclar Vivadent)

Self-cure & dual-cure luting composite, with light curing option

Kit

\$199.95

9g syringe, Monobond Plus bottle, 15 mixing tips, 50 applicator brushes

Transparent 246-4060

MO0 246-4061

MO1 246-4062

Refills

\$129.95

1 x 9g syringe

Transparent 246-4063

MO0 246-4064

MO1 246-4065



Premier® Implant Cement (Premier)

Ideal for both implant retained crown & long-term temporaries

Kit

5ml automix syringe, & 10 tips

246-4070 **\$49.95**



Luting Cement (ibt Med)

Automix dual-cure

\$29.50

Kit

10g Syringe

(0120-7072) A2 242-5040

(0120-7073) A4 242-5041

(0120-7074) Blue 242-5042

(0120-7075) White Opaque 242-5043



Fax your order to: 310.537.1502



Permanent Cement

veneer

Accolade™ PV (Danville)

Light-cured thixotropic flowable composite, high viscosity

Veneer Kit 282-2120 **\$299.95**
5g Translucent, 3g Light, 3g Extra Light, 3g White Opaquer, 3g Brown & 3g Yellow (in PV Composite), 5g Translucent, 3g light, 3g Extra Light, 3g White Opaquer, 3g Brown & 3g Yellow (in Try-In, ultra slow set composite), Prelude SE Bonding system, SereEtch Gel, S-Bond



	PV Composite Refills	Try-In Refills
	\$26.95	\$26.95
Light 3g	282-2121	282-2131
Extra Light 3g	282-2122	282-2132
White Opaquer 3g	282-2123	282-2133
Brown 3g	282-2124	282-2134
Yellow 3g	282-2125	282-2135
Translucent, 5g: \$44.95	282-2126	282-2136

Accessories

Prelude SE Self Etch Kit 282-2260 **\$119.95**

SureEtch™ **\$19.95**

Bottle 5ml 275-2262

Syringe 5ml 275-2263

S-Bond 275-2052 **\$34.95**

RelyX™ Veneer Cement (3M ESPE)

Veneer cement provides simple & efficient bonding of porcelain, composite, & ceramic veneers.

Trial Kit 7613TK 242-3620 **\$199.95**

3g syringe Translucent Veneer Cement, 2g syringe Try-in paste (Translucent), 5ml ceramix Primer, 3ml etchant, 25 Etchant dispensing tips, 20 applicators, technique card

Kit 7613 242-3621 **\$449.95**

5 x 3g syringes of Veneer Cement & 5 x 2g, Try-In Paste Syringes, 4.5g ARC Adhesive Resin Cement Paste Shade Universal in Clicker Dispenser, 6 ml Adper Single Bond Dental Adhesive, 5 ml Ceramic Primer, 3 ml Scotchbond Phosphoric Etchant Syringe, & accessories

Refills

Syringes, 3g **\$59.95**

Try-in Paste, 3g **\$49.95**

	Syringe	Try-in Paste
Translucent	242-3650	242-3660
White Opaque	242-3651	242-3661
A1/Light Yellow	242-3652	242-3662
A3/Opaque Yellow	242-3653	242-3663
A5/Dark	242-3654	242-3664
B0.5/White	242-3655	242-3665



Mojo™ Veneer Cement (Pentron)

The most biocompatible material for soft tissue & teeth. Dual-cure hydrophylic bis-GMA

Kit 242-5370 **\$199.95**

2.3g cement, 16g Try-in paste, 4ml Bond-1 Primer/Adhesive, 1ml syringe 37% Phosphoric acid etching gel, Silane, 40 Elephant tips

Refills

Cement - 2.3g syringe, 5 tips, instructions **\$29.95**

Try-In Gel - 1.6g syringe, 5 tips, instructions **\$14.95**

	Cement	Try-In Gel
Light	242-5371	242-5375
Dark	242-5372	242-5376
Clear	242-5373	242-5377
Movie Star White	242-5374	242-5378



NX3 (Kerr)

Third generation, truly universal perfection in permanent cements

Intro Kit 33642B

Dual cure 5g syringe of each shade: clear, white and yellow; 1.8g syringe light-cure of each shade, 1 Syringe Try-in Gel 3g of each shade, 1 Bottle Silane Primer 5ml; 30 Automix Tips; technique guide

242-4100 **\$349.95**

Refill

Automix Dual-Cure Syringe, 5g

Try-In Gel 3g

	Automix	Try-In Gel
	\$99.95	\$33.75
Clear	242-4101	242-4106
White	242-4102	242-4107
Yellow	242-4103	242-4108
Bleach	242-4104	242-4109
White Op	242-4105	242-4110



Kleer-Veneer™ (Pulpdent)

Light-cure veneer cement

Kit **\$72.95**

4 x 1.2ml syringe with 20 tips 242-4410

Syringe **\$24.95**

	Syringe
1.2ml with 10 tips	
Clear	242-4411
Cream	242-4412
Pink	242-4413
White	242-4414



Permanent Cement

primer

RelyX™ Ceramic Primer (3M ESPE)

Silane-based primer, 5ml

242-3657 **\$99.95**



Clearfil™ Ceramic Primer (Kuraray)

One-step silane coupling agent

4 ml Ceramic primer 242-3705 **\$69.95**



Clearfil™ Porcelain Bond Activator (Kuraray)

Silane coupling agent

Activator, 4 ml 242-3740 **\$59.95**



Alloy Primer (Kuraray)

Used to secure the adhesive capability of resins to gold

Alloy Primer, 4 ml 242-3745 **\$54.95**



Copaliner (Bosworth)

Dentin resin varnish & sealant.

Complete Kit 246-1010 **\$12.95**

14 ml each varnish & thinner

Varnish, 14ml 246-1011 **\$7.95**

Thinner, 14ml 246-1012 **\$5.95**



Cavity Liners

self-cure

Dycal (Caulk)

Cavity liner/base, radiopaque calcium hydroxide composition

Standard Pack **\$39.95**

13g base + 11g catalyst + mixing pads

Dentin 246-1700

Ivory 246-1720



Pulp Capping Paste (Pulpdent)

The original pre-mixed calcium hydroxide dressing, radiopaque

Kit 246-1750 **\$29.95**

3ml syringe + 24 sterile dispensing needles



Life® (Kerr)

Base & catalyst package, self-cure

Regular Set 246-1770 **\$49.95**

Life Fast Set 246-1775 **\$49.95**



Cavity Liners

liquid

Copalite (Cooley & Cooley)

Kit 246-1005 **\$17.95**

Solvent & Thinner, 0.5 oz

246-1000 **\$11.95**

Varnish, 0.5 oz 246-1006 **\$14.95**



Copal Cavity Varnish (Deepak)

0.5oz

Standard 246-1050 **\$3.50**

w/Fluoride 246-1051 **\$3.95**




Hydrox (Bosworth)

Standard Package **\$29.50**

15g base & 13g catalyst, mixing pad

Dentin 246-1780

Ivory 246-1781



CEMENT & LINERS

Ketac-Bond (3M)

Radiopaque glass ionomer base & liner

Powder Yellow, 10g	246-1800	\$59.95
Liquid, 10ml	246-1801	\$39.95



Ketac-Bond Aplicap (ESPE)

Capsulated radiopaque glass ionomer base

Yellow Caps (50)	246-1810	\$199.50
Aplicap Activator	246-1811	\$49.50



Bio-Cap™ (DenMat)

Adhesive base-liner & pulp cap, fast, easy & efficient method of pulp capping with a formula that has clinical biocompatibility. Self-cure

Perio/Endo Kit A3.5	246-1850	\$119.95
Syringe, 4ml	246-1855	\$54.95



Fuji Lining LC Paste (GC America)

1 Cartridge refill, dispenser, accessories

Starter pkg	246-1930	\$129.50
--------------------	----------	-----------------

1 Paste pack cartridge, mixing pad, dispenser & spatula



Paste pkg Refill	246-1931	\$99.50
-------------------------	----------	----------------

7g Paste pack cartridge, & mixing pad

Paste Pack Dispenser	240-3705	\$59.95
-----------------------------	----------	----------------

Vitrebond (3M ESPE)

Light-cured glass ionomer liner/base material

Introductory Kit	242-3250	\$199.95
-------------------------	----------	-----------------

9g powder & 5.5ml liquid

Liquid, 5.5ml	242-3251	\$129.95
Powder, 9g	242-3252	\$99.95



Vitrebond Plus (3M ESPE)

Light-cured glass ionomer liner/base packs

Single pkg	242-3260	\$149.95
-------------------	----------	-----------------

Includes clicker dispenser, large mixing pad, instructions

Double pkg	240-3265	\$259.95
-------------------	----------	-----------------

Includes clicker dispenser, large mixing pad, instructions



Cavity Liners

light-cure

Lime-Lite (Pulpdent)

Kit	246-1900	\$29.50
------------	----------	----------------

4 x 1.2ml syringes + 8 tips



Biner LC (Meta Biomed)

Hydroxy calcium phosphate, with fluoride

Kit	246-1905	\$29.50
------------	----------	----------------

2 x 2ml syringes + 10 tips



Fuji Bond LC (GC America)

Light-cured resin reinforced glass ionomer dentin/enamel bonding-liner.

Standard Pack	246-1920	\$199.50
----------------------	----------	-----------------

5g powder, 7 ml liquid, 5.7 ml GC cavity conditioner & accessories

Powder, 5g	246-3075	\$129.50
Liquid 7ml	246-3076	\$89.50



Fuji Lining LC (GC America)

Light-cured glass ionomer cement for lining. Radiopaque & fluoride releasing

1:1 Package	246-2000	\$139.50
--------------------	----------	-----------------

Contains 10g powder, 6.8 ml liquid, & measuring scoop

Powder, 10g	246-2010	\$79.50
Liquid, 6.8ml	246-2011	\$69.50



Desensitizers

Gluma Desensitizer (Heraeus)

5 ml bottle	246-2200	\$109.95
Clinic Pack - 3 x 5 ml bottle	246-2210	\$279.95
Single Dose - 40 x 0.075ml	246-2215	\$69.95



Gluma Desensitizer PowerGel

(Heraeus)

Kit

4 x 1gm syringe, 20 brush cannulas and pictorial card 246-2500 **\$109.95**

Brush Refills

60/pkg 246-2501 **\$14.95**



Gluma Comfort Bond & Desensitizer

(Heraeus)

5th generation adhesive, 4-META strength, exceptional bond strength, adhesion to metal and dentin. Suitable for both moist & dry bonding. Direct filling light curing composite

Single Pack, 4ml 246-2220 **\$109.95**

Economy Pack, 3 x 4ml 246-2221 **\$299.95**

Single Dose, 40 x 0.1ml 246-2222 **\$99.95**



Pac Desensitizer (Pac-Dent)

1-step solution to minimize dentin or post operative sensitivities

10ml bottle 246-2300 **\$27.50**



MicroPrime™ (Danville)

Universal desensitizing agent

MicroPrime B 246-2310 **\$49.95**

35% HEMA, fluoride & benzathonium
10ml bottle

MicroPrime G 246-2315 **\$49.95**

35% HEMA, and 5% glutaraldehyde
10ml bottle



Dentin Desensitizer (Pulpdent)

12 ml bottle 246-2350 **\$34.99**



Hurri-Seal (Beutlich)

Antibacterial water-based solution

12 ml bottle 246-2400 **\$39.50**



Bracket Adhesive

self cure

Fuji Ortho P&L (GC America)

Self-Cured Ortho Bonding Adhesive
Powder & liquid

Standard Kit 439503

40g powder, 2 x 6.7ml bottles of liquid, scoop, mixing pad 246-2910 **\$199.95**

Liquid (6.7ml) 209 246-2915 **\$49.95**

Powder (40g) 210 246-2916 **\$129.95**



Fuji Ortho Band Paste Pack (GC America)

Blue-colored resin modified glass ionomer for orthodontic bands, automix

Starter Kit 439460 246-2930 **\$129.50**

1 refill box, mixing pad, dispenser, & spatula

Refill Kit 001888 246-2935 **\$79.50**

2 x 13g cartridges

Paste Pack Dispenser 240-3705 **\$59.95**



Glass Ionomer Band Cement (Pulpdent)

Provides strong bond to metal bands, radiopaque, fast hardening

Kit 246-2940 **\$39.95**

30g powder, 15ml liquid, scoop, mixing pad

Bulk Pack 246-2945 **\$169.95**

240g powder, 120ml liquid, scoop, mixing pad

Powder

30g 246-2952 **\$34.95**

240g 246-2953 **\$129.95**

Liquid

15ml 246-2954 **\$7.50**

120ml 246-2955 **\$59.95**



(Dentonics/Master-dent)

Bracket Adhesives

Light Cure Adhesive

4 x 5g syringe, 7g primer, 15ml etch, spatulas, pads, brushes, brush holder

246-2960 **\$59.95**



2-Paste Self Cure Adhesive

18g catalyst paste, 18g base paste, 7g catalyst resin, 7g base resin, 9g etch, spatulas, mixing pads

246-2965 **\$49.95**



Bracket Adhesive

light/dual cure

Fuji Ortho Band LC (GC America)

A blue resin modified glass ionomer for bonding orthodontic bands, can also be used for the cementation of appliances retained with stainless steel crowns



Starter Kit 439466

Automix, 3.5g Paste Pak Cartridge, 2 -1.4 g ortho gel conditioner syringe, 20-mixing tips, 15 x brush tips 246-3000 **\$199.95**

Refill - Paste Pak 003105

2 x 13.5 g Paste Pak cartridge 246-3010 **\$155.95**

Fuji Ortho LC - P&L (GC America)

Powder & Liquid, A light-cured reinforced glass ionomer cement

Introductory Kit 439401

15g powder, 8ml liquid & Accessories 246-3050 **\$119.95**



Standard kit 439403

40g powder, 2 x 8ml liquid & Accessories 246-3055 **\$249.50**

Powder (40g) 000029 246-3056 **\$149.95**

Liquid (8ml) 000030 246-3057 **\$59.95**

Fuji Ortho LC - Caps (GC America)

Capsule Package

Starter Kit 439451

50 capsules, applier, & Ortho conditioner 246-3060 **\$299.50**



Refill Package #439450 50 capsules 246-3065 **\$229.50**

Capsules Applier IV 246-3066 **\$69.95**

Fuji Ortho LC - Automix

Automix Package

Automix Starter Kit 439487

Cartridge 13.3g, Paste Pak Dispenser, 2 Conditioner Syringes 1.4g each, accessories 246-3070 **\$299.95**



Refill Package 003103 2 Paste Pak Cartridges

(13.3g) 246-3075 **\$229.95**

Band-Rite™ (Pulpdent)

Fluoride releasing, dual-cure paste

Standard Kit

2 x 12g syringe base, 2 x 12g syringe catalyst, accessories 246-3080 **\$89.95**

Refill

12g syringe base, 12g syringe catalyst 246-3085 **\$39.95**



Heliosit® Orthodontic (Ivoclar Vivadent)

Bracket bonding material, light-cured bonding material for the direct bonding of orthodontic brackets.

- Use with plastic or metal brackets.
- Low film thickness
- High bond to enamel
- Three year shelf life

3 x 2.5g syringes 246-3120 **\$89.95**



Aegis® (Bosworth)

The worlds first light-cured orthodontic adhesive containing amorphous calcium phosphate (ACP) It is used for bonding orthodontic brackets and bands to etched enamel.

Standard Package

3 x 2.5g syringes 246-3090 **\$169.95**



OBA™ (Pulpdent)

Light-cure fluoride relasing adhesive

2-Syringe Kit 246-3100 **\$49.95**

2 x 4g syringes, 5ml resin bonding agent, 6ml Etch-Rite gel, accessories

OBA Syringe only, 4g 246-3105 **\$14.95**

Bonding resin, 5ml 246-3106 **\$14.95**



Band-Rite™ (Pulpdent)

Fluoride releasing, dual-cure paste

Standard Kit

2 x 12g syringe base, 2 x 12g syringe catalyst, accessories 246-3080 **\$89.95**

Refill

12g syringe base, 12g syringe catalyst 246-3085 **\$39.95**



OBA™ Pre-filled (Pulpdent)

20 pre-filled tips of bonding agent, 6ml Etch-Rite gel, accessories

246-3110 **\$59.95**

