

EVACUATION PRODUCTS

Evacuation Tips - HVE

plastic

(Unipack)

Evacuator Tips		\$3.95
Combo vented/non-vented 100/Bag		
White	430-7585	
Blue	430-7586	



(StarryShine)

Vented/non-vented		
1000/case	430-7520	\$34.95



(Crosstex)

Combo-Tip®		
Soft tip evacuator vented on one end		
(Same as Hygovac)		
100/bag	430-7530	\$3.95



Maxvac Plus Combo Tip®

100/bag		\$4.95
Blue	430-7551	Yellow 430-7554
Green	430-7552	Purple 430-7555
Pink	430-7553	
Mix - case - 2 bags of		
each color	430-7556	\$45.95



(Plasdent)

Spectrum		
Bendable, 7 assorted colors		
100/bag	430-7535	\$7.95



Pastel

Bendable, 6 assorted colors		
100/bag	430-7540	\$7.95



Hi-Vac Multi Colors

100/bag	430-7545	\$3.95
---------	----------	--------



(Plasdent)

Oral Evacuator



8008 8015 8018 8888 8048 8048A 8028

Bendable-compare to Monoject
pkg of 100

8008 - Hi Vac	430-7561	\$3.95
8015 - OraJet	430-7562	\$3.95
8018 - OraVac	430-7563	\$4.95
8888 - DualVac	430-7564	\$3.95
8048 - TrkOVac	430-7565	\$7.95

Autoclavable

8028 - Ultra Vent, 50/Bag	430-7566	\$6.95
8048A - TrkOVac, 25/Bag	430-7567	\$7.95

Tip Inserts

pkg of 25		\$6.50
8015 - TP, tip inserts	430-7571	
8015 - HS, surge inserts	430-7572	
8015 - HJ, jet inserts	430-7573	



TP HS HJ

Oratip III (Pascal)

Fits all high volume evacuation systems		
50/Bag	430-7550	\$30.95



Accessories

Combo-Tip® Bender (Crosstex)

430-7600	\$8.95
----------	--------



HVE Tips Screen

100/pk		
Vista - Blue	430-7601	\$12.95
Plasdent - Clear	430-7602	\$6.95



EVACUATION PRODUCTS

Evacuation Tips - HVE

metal

(Vista)

High Volume Aspirators

Stainless steel, autoclavable

Slotted Tip (Reduces Tissue-Grab) \$28.95

802001: Fits HVE Cut-off Valve
Equivalent to 46RSA 435-2001

802002: 41/4" Total Length
Equivalent to 44RSA 435-2002

802005: Flattened Surface
6 x 4mm opening
Equivalent to 46FSA 435-2005

802006: Flattened Surface
8 x 5mm opening
Equivalent to 56FSS 435-2006

802007: Flattened Surface
9.5 x 6.5mm opening
Equivalent to 66FSS 435-2007

802008: Flattened Surface
11.5 x 7.5mm opening
Equivalent to 76FSSS 435-2008

Back-Cut (Reverse Bevel) \$28.95

802009: 6 x 4mm opening
Equivalent to 46FSAABC 435-2009

802010: 8 x 5mm opening
Equivalent to 56FSSABC 435-2010

802011: 9.5 x 6.5mm opening
Equivalent to 66FSSBC 435-2011

802012: 11.5 x 7.5mm opening
Equivalent to 76FSSSBC 435-2012

Short (41/4" Long) \$28.95

802013: 6 x 4mm opening
Equivalent to 46FSA 435-2013

802014: 8 x 5mm opening
Equivalent to 54FSS 435-2014

802015: 9.5 x 6.5mm opening
Equivalent to 64FSS 435-2015

802016: 11.5 x 7.5mm opening
Equivalent to 74FSSS 435-2016

Universal (10 x 13.5mm)

802003: Fits HVE Cut-off Valve \$28.95
10 x 13.5mm opening
Equivalent to 46STR 435-2003

802004: Fits HVE Cut-off Valve \$28.95
10 x 13.5mm opening
Equivalent to 76R 435-2004



Straw Aspirator

802017: Features Side Venting \$74.95

Suction Adjustment Holes

Equivalent to Tip-A-Dilly 437-2100

Straws, 200/pkg 437-2101 \$10.95



Accessories

Suction Tube Adaptors (Pac-Dent)

11mm to 6.5mm 437-0010

11mm to 11mm 437-0011

11mm to 16mm 437-0012



HVE Adapters (Plasdent)

1/4" surgical tubing

Blue 438-1600 Red 438-1603

Yellow 438-1601 Black 438-1604

Green 438-1602



Saliva Ejectors

(Unipack)

Disposable Saliva Ejectors

Soft non-removable tip, vented tip

100/pkg \$3.95

White 436-1100

Clear, blue tip 436-1101

Pink, pink tip 436-1102

Blue, clear tip 436-1103

Green, clear tip 436-1104



(Starry Shine)

1000/case, 5 3/4" \$34.95

White 436-1000

Clear, blue tips 436-1001



(Plasdent)

100/case, 5 3/4" \$3.95

Flexjet™

White, white tip 436-1010

Clear, white tip 436-1011

white, blue tip 436-1012

Clear, blue tips 436-1013



EVACUATION PRODUCTS

Saliva Ejectors

Flexjet™ \$3.95
 100/case, 5 3/4"
 Clear multicolor 436-1017
 Opaque multicolor 436-1018



(Plasdent)

Flexform
 100/case, 5 3/4"
 White w/clear adapter 436-0030 \$4.95



(Crosstex)

Non-Removable Tip
 Comfort Plus™ Premium \$4.95
 100/bag
 White/White 436-0040
 Clear/White 436-0039



Comfort Plus™ Premium - Scented
 100/bag
 Bubble Gum - White/Pink 436-0037 \$4.95
 Mint - White/Green 436-0038 \$4.95



Premium®
 100/bag, 5 3/4" \$4.95
 Clear/Blue 436-0020
 Clear/White 436-0021
 White/White 436-0022
 Pink/Pink 436-0023
 Multicolor 436-0024



Advantage
 100/bag, 5 3/4" \$3.49
 Clear/Blue 436-0025
 Clear/Clear 436-0026
 White/White 436-0027
 Blue/Blue 436-0028

Free Form / Hygoforming

Saliva ejector/tongue retractor
 White, 100/bag 436-0035 \$9.95
 Adapters, 10/bag 436-0036 \$9.95



SE Cushion (Zirc)

Used in conjunction with the saliva ejector
 100/Pack 436-1200 \$18.95



Ligua-Fix (Zirc)

Holds tongue in a steady position, helps protect the tongue from injuries
 10-Pack 436-1150 \$12.95
 48-Pack 436-1151 \$44.95



Saliva Ejectors

metal

(Vista)

30° Angle \$25.95
 Tapered Tip (Reduces clogging)
 806001: 1.5mm opening 436-1250
 Equivalent to 14P1-1/4
 806002: 2.5mm opening 436-1251
 Equivalent to 14P2-1/4
 806003: 3.0mm opening 436-1252
 Equivalent to 14P3-1/4
 Equivalent to 15P3-1/4



90° Angle \$25.95
 Tapered Tip (Reduces clogging)
 806004: 1.5mm opening 436-1253
 Equivalent to 14P1-90-1/4
 806005: 2.5mm opening 436-1254
 Equivalent to 14P2-90-1/4
 806006: 3.0mm opening 436-1255
 Equivalent to 14P3-90-1/4



Surgical Tips

plastic

(Unipack)

Disposable Surgical Aspirator Tips
 25/bag \$4.95
 X-Small Blue 1/16" 435-0500
 Small White 1/8" 435-0501
 Large Green 1/4" 435-0502



EVACUATION PRODUCTS

(StarryShine)

Disposable, pkg of 25	\$4.95
X-Small - Blue 1/16"	435-0005
Small - White 1/8"	435-0006
Large - Green 1/4"	435-0007
X-large - Blue 5/6"	435-0008



(Pac-Dent)

Disposable or Autoclavable		
25/pkg		\$6.50
Green, large	435-0010	
White, small	435-0011	



(Vista)

25/pkg		\$6.95
Small, White 1/8"	435-1010	
Medium, Green 1/4"	435-1011	
Large, Blue 3/8"	435-1012	



(Plasdent)

Oralsurge II		\$6.95
25/Bag		
1/4" Orifice		
Aquamarine	435-0021	
White	435-0022	
1/8" Orifice	435-0023	
White		



Oralsurge III		\$6.95
3/8" Orifice	435-0024	
vented tip, Seablue		



Oralsurge-ST		\$7.95
1/16" Orifice		
Blue	435-0025	

Oratip II (Pascal)		\$14.95
50/pkg	435-0026	

Chromasurge™		\$11.95
1/8" orifice, 100/Bag		
Assorted	435-0027	



Surgical Aspirators

metal

(Vista)

Stainless steel, autoclavable
Fits HVE Cut-Off Valve



3/16" Dia Tubing

Tapered Tip (Reduces clogging) \$28.95

801001: 1.5mm opening
Equivalent to 15P1A 437-0500

801002: 2.5mm opening
Equivalent to 15P2A 437-0501

801003: 3.0mm opening
Equivalent to 15P3A 437-0502

Slotted Tip (Reduces tissue-grab) \$28.95

801004: 4.0mm opening
Equivalent to 16RSA 437-0505

Right Angled Tapered Tip (Reduces clogging) \$28.95

801005: 1.5mm opening
Equivalent to 15P1-90A 437-0510

801006: 2.5mm opening
Equivalent to 15P2-90A 437-0511

801007: 3.0mm opening
Equivalent to 15P3-90A 437-0512

Right Angled Slotted Tip (Reduces tissue-grab)

801008: 4.0mm opening
Equivalent to 16RS-90A 437-0515 \$28.95

1/4" Dia Tubing \$28.95

Tapered Tip (Reduces clogging)

801009: 1.5mm opening
Equivalent to 46P1A 437-0520

801010: 2.5mm opening
Equivalent to 46P2A 437-0521

801011: 3.0mm opening
Equivalent to 46P3A 437-0522

5/6" Dia Tubing \$28.95

Tapered Tip (Reduces clogging)

801012: 2.5mm opening
Equivalent to 56P2 437-0530

801013: 3.0mm opening
Equivalent to 56P3 437-0531

801014: 4.0mm opening
Equivalent to 56P4 437-0532

Frazier Style \$42.95

Fits 1/4" Surgical Tubing

805001: 1.5mm opening
Equivalent to FR1 437-0535

805002: 2.5mm opening
Equivalent to FR2 437-0536

805003: 3.0mm opening
Equivalent to FR3 437-0537



EVACUATION PRODUCTS

Evacuation System Cleaners

liquid

E-Max (FCI)
Enzymatic super concentrated
1 Gallon 434-1030 \$19.95



Pro E-Vac (Certol)
Concentrated
1 Gallon w/pump 434-1040 \$39.95
1oz Unit dose (24/box) 434-1041 \$16.95



SaniVac (Crosstex)
Concentrated
1 Gallon w/pump 434-1042 \$54.95



Vacu Blast (Premier)
Enzymatic Concentrated
32 fl. oz. 434-1043 \$54.95



Vacuosol Ultra (Biotrol)
Starter Kit 434-1044 \$108.95
Refill
32 oz Bottle 434-1045 \$59.95
96 oz Bottle 434-1046 \$121.95
1-2-3 Atomizer 434-1047 \$52.95



NeutraVac (Biotrol)
32 oz Bottle 434-1051 \$61.95
96 oz Bottle 434-1052 \$125.95



PureVac (Sultan Chemists)
2 Liter Bottle 434-0055 \$61.95
5 Liter Bottle 434-0056 \$129.95
Walkabout Measuring & Dispensing Unit 434-0059 \$87.95



PureVac SC (Sultan Chemists)
1 oz makes 1 quart of solution
1 oz pump, non-Foaming & Biodegradable



21132 2-Liter (67oz) Bottle (67 treatments) 434-0065 \$63.95
21135 5-Liter (169 oz) Bottle (169 treatments) 434-0066 \$139.95


Evacuation System Cleaners

powder

E.S.C. Powder (California Dental)
2 lb powder 434-1000 \$12.95



Citrizyme™ (Pascal)
Dual enzymatic cleaner for ultrasonic and evacuation systems, ultra concentrated
300g powder 434-1020 \$34.95
Unit Dose, 50/pkg 434-1025 \$37.95
900g powder 434-1026 \$76.95



VAC Attack (Premier)
800g powder 434-1070 \$52.95



Evacuation Cups

Liners for A-DEC
6oz white
250/pkg 433-1000 \$19.95
Case 10 x 250/pkg 433-1020 \$140.95



Dry Oral Cup (Plasdent) \$3.95
Beige, autoclavable, single pack
Type I - 8118 433-1010
Type II - 8119 433-1015



EVACUATION PRODUCTS

Disposable Traps


plastic



8218 8178 8114 8158 8334

Screen-4-Vac (Plasdent) Clear/White \$20.95
144 Traps/box Colored \$23.95

8218, Clear 438-0100
B8218 - Blue 438-0101
Fits A-dec w/plastic canister, 2 1/8" Dia



8178, Clear 438-0102
G8178, Green 438-0103
fits Belmont, Biotec, Chayes, Den-Tal-Ez, Dentech,
Pelton & Crane, Proma, 1 7/8" Dia, Clear

8114, White 438-0104
Universal, fits most cuspidor bowls

8158, White 438-0105
fits A-dec Casade, Performer cuspidors bowls, 1 5/8" Dia

8334, Clear 438-0106
fits Pelton & Crane Spirit with white plastic canisters, 3 3/4" Dia



8814 8214 8078 8112 8100

Clear/White \$20.95 Color \$23.95

8814 Clear 438-0017
fits A-dec Casade & Performer, 1 1/4" Dia



8214 Clear 438-0018
Y8214 Yellow 438-0019
fits A-dec w/metal canister, 2 1/4" Dia

8078 White 438-0020
fits Siemens, Planmeca & many cuspidors with small openings
Small

8112 Clear 438-0021
fits DCI/Belmont Reality models, System 6000 models, 1 1/2" Dia

8100 438-0022 \$9.95
Saliva Ejector traps, 72 pcs/bag, one size fits all

Evacuation Screens (Palmero)

Disposable Traps & Strainers

pkg of 144 \$42.95

1805 Blue, Dia, 2 1/8" 438-0001
Fits: A-dec with
plastic canister



1806 Yellow, Dia, 2 1/4" 438-0002
Fits: A-dec with
metal canister



1807 Green, Dia, 1 7/8" 438-0003
Fits: Belmont, Biotech,
Pelton & Crane, Proma



Dispos-A-Screen Trap (Zirc)

Saliva Ejector Screens

25/pack \$15.95

144/pack \$61.95



25-pack 144-pack

Mauve 1-1/4" 438-1500 438-1510

Yellow 2-1/4" 438-1501 438-1511

Blue 2-1/8" 438-1502 438-1512

Coral 2-5/8" 438-1503 438-1513

Beige 2-7/8" 438-1504 438-1514

Gray 2-3/4" 438-1505 438-1515

Saliva Ejector Screens (Zirc)

25/pack 438-1518 \$9.95



Cuspidor Screens (Zirc)

25/pack \$17.95

144/pack \$69.95

Small 1-5/8"

Fits Adec, Cascade and Performer cuspidor
bowls, gray



Large universal 2-1/2"

Fits Adec, Ampco, Belmont, Dab Atante,
Dansereau, DCI Den-Tal-Ez, Forest, Knight,
Marus, MDT, Proma, Ritter, Schein



25/pack Small Large 438-1520 438-1525

144/pack 438-1521 438-1526



EVACUATION PRODUCTS

Disposable Traps

plastic

ProEvac™ Disposable Traps (Pac-Dent)
144 filters/pack

\$32.95



PD-5500 Dia 1 1/4"
fits A-Dec, Cascade & Performer

438-1700

5500

PD-5501 Dia 2 1/8"
fits A-Dec w/plastic canisters

438-1701



5501

PD-5502 Dia 2 1/4"
fits A-Dec w/metal canisters

438-1702



5502

PD-5503 Dia 1 7/8"
fits Belmont, Biotec, Chayes,
Den-Tal-Ez, Dentech, Pelton & Crane,
Proma, etc

438-1703



5503

PD-5505 Dia 3 3/4"
fits Pelton & Crane Spirit w/white
plastic canisters

438-1704



5505

PD-5508 Dia Universal
fits all saliva ejectors, 72/box

438-1705



5508

PD-5509 fits DCI/Belmont
Realisty models, System 6000 models

438-1706



5509

PD-5510 Dia 3"
fits Sirona - C2, C3, C4, C6, C8,
Pelton & Crane S1, S2

438-1707



5510

PD-5512 Dia 2 3/4"
fits A-dec 500

438-1708



5512

PD-6100 Dia Universal
fits most Cuspidor Bowls
including: Adec, Ampco, Belmont,
Dabi Atlante, Dansereau, DCI,
Den-Tal-Ez, Forest, Knight, Marus,
MDT, Proma, Ritter, Schein

438-1709



6100

PD-6200 Dia Universal
fits Siemens, Planmeca, and
many cuspidors with small
diameter openings

438-1710



6200

PD-6500 Dia Universal
fits Kavo the environment

438-1711



6500

Disposable Canisters

holder



Disposable Canister (Plasdent)

Each

\$9.95

DC8-2400, fits Dentsply-MVS 0.75, 1.5, 4.0,
before 1991 (4 1/2" W x 5 1/8" H)

438-0030

DC8-2300, fits Den-Tal-Ez dual after 1978,
Dentsply after 1991, Matrix 1000, 2000 &
Dynamix models (3 1/2" W x 4 3/8" H)

438-0031

DC8-2350, fits Apollo 2 hp, 3 hp single &
tandem, Air Techniques 1997 & later dual
models, Vacudent-Workhorse dual & newer
Vacstar (3 1/2" W x 4 3/8" H)

438-0032

DC8-2200, fits Air Techniques-Vacstar, Apollo
1 hp 1.5 single, Tandem, 1, 2 hp Den-Tal-Ez
single, Ohmeda-matrix, Turbine Industries,
Techwest-Whirl Vac, Vacudent-Workhorse
single, plus & other units. (2 3/4" W x 3 5/8" H)

438-0033

Easy-E-Trap (Pac-Dent)

Kit consist of waste canister & lid, filter
bowl, filter screen, replacement rubber gasket,
& warning sticker



Single

Box

\$7.95

\$43.95

Model #2200 438-2000

438-2005 12/box

Model #2300 438-2001

438-2006 8/box

Model #2350 438-2002

438-2007 8/box

Model #2400 438-2003

438-2008 8/box

Model #2500 438-2004

438-2009 8/box

Mesh Filters (Pac-Dent)

Filters for central vacuum pump
1/pkg



H: 3.75", dia 1.50"
fits model #2200

438-0050

\$2.95

H: 4.25", dia 2.50"
fits model #2300/2350

438-0051

\$3.49

