

# EVACUATION PRODUCTS

## Evacuation Tips - HVE

plastic

(Unipack)

### Evacuator Tips

Combo vented/non-vented

100/Bag 430-7585 **\$2.95**



(StarryShine)

Vented/non-vented

1000/case 430-7520 **\$24.95**



(Winner)

Vented/non-vented

1000/case 430-7525 **\$24.95**

(Crosstex)

### Combo-Tip®

Soft tip evacuator vented on one end  
(Same as Hygovac)

100/bag 430-7530 **\$3.95**



### Maxvac Plus Combo Tip®

100/bag **\$4.49**

Blue 430-7551 Yellow 430-7554

Green 430-7552 Purple 430-7555

Pink 430-7553

Mix - case - 2 bags of  
each color 430-7556 **\$39.95**



(Plasdent)

### Spectrum

Bendable, 7 assorted colors


100/bag 430-7535 **\$6.95**



### Pastel

Bendable, 6 assorted colors

100/bag 430-7540 **\$6.95**



### Hi-Vac Multi Colors

100/bag 430-7545 **\$3.75**



## Oral Evacuator



8008 8015 8018 8888 8048 8048A 8028

Bendable-compare to Monoject  
pkg of 100

8008 - Hi Vac 430-7561 **\$2.75**

8015 - OraJet 430-7562 **\$2.75**

8018 - OraVac 430-7563 **\$3.50**

8888 - DualVac 430-7564 **\$3.50**

8048 - TrkOVac 430-7565 **\$7.95**

### Autoclavable

8028 - Ultra Vent, 50/Bag 430-7566 **\$6.95**

8048A - TrkOVac, 25/Bag 430-7567 **\$7.95**

### Tip Inserts

pkg of 25 **\$6.50**

8015 - TP, tip inserts 430-7571

8015 - HS, surge inserts 430-7572

8015 - HJ, jet inserts 430-7573



TP HS HJ

### Oratip III (Pascal)

Fits all high volume evacuation systems

50/Bag 430-7550 **\$24.95**



## Accessories

### Combo-Tip® Bender (Crosstex)

430-7600 **\$4.50**



### HVE Tips Screen

100/pk

Vista - Blue 430-7601 **\$12.95**

Plasdent - Clear 430-7602 **\$6.95**



### Suction Tip Sleeves (Pac-Dent)

Fully autoclavable. Blue color

6 tips/bag 437-0020 **\$9.95**



## Evacuation Tips - HVE

metal

(Vista)

### High Volume Aspirators

Stainless steel, autoclavable

#### Slotted Tip (Reduces Tissue-Grab) **\$21.95**

**802001:** Fits HVE Cut-off Valve  
Equivalent to 46RSA 435-2001

**802002:** 4 1/4" Total Length  
Equivalent to 44RSA 435-2002

**802005:** Flattened Surface  
6 x 4mm opening  
Equivalent to 46FSA 435-2005

**802006:** Flattened Surface  
8 x 5mm opening  
Equivalent to 56FSS 435-2006

**802007:** Flattened Surface  
9.5 x 6.5mm opening  
Equivalent to 66FSS 435-2007

**802008:** Flattened Surface  
11.5 x 7.5mm opening  
Equivalent to 76FSS 435-2008



#### Back-Cut (Reverse Bevel) **\$21.95**

**802009:** 6 x 4mm opening  
Equivalent to 46FSABC 435-2009

**802010:** 8 x 5mm opening  
Equivalent to 56FSSABC 435-2010

**802011:** 9.5 x 6.5mm opening  
Equivalent to 66FSSBC 435-2011

**802012:** 11.5 x 7.5mm opening  
Equivalent to 76FSSSBC 435-2012



#### Short (4 1/4" Long) **\$21.95**

**802013:** 6 x 4mm opening  
Equivalent to 46FSA 435-2013

**802014:** 8 x 5mm opening  
Equivalent to 54FSS 435-2014

**802015:** 9.5 x 6.5mm opening  
Equivalent to 64FSS 435-2015

**802016:** 11.5 x 7.5mm opening  
Equivalent to 74FSS 435-2016



#### Universal (10 x 13.5mm)

**802003:** Fits HVE Cut-off Valve **\$19.95**  
10 x 13.5mm opening  
Equivalent to 46STR 435-2003

**802004:** Fits HVE Cut-off Valve **\$19.95**  
10 x 13.5mm opening  
Equivalent to 76R 435-2004



### Straw Aspirator

**802017:** Features Side Venting **\$39.95**  
Suction Adjustment Holes  
Equivalent to Tip-A-Dilly 437-2100  
Straws, 200/pkg 437-2101 **\$7.95**



### Accessories

#### Suction Tube Adaptors (Pac-Dent)

11mm to 6.5mm 437-0010 **\$9.50**  
11mm to 11mm 437-0011 **\$11.50**  
11mm to 16mm 437-0012 **\$12.50**



#### HVE Adapters (Plasdent) **\$2.95**

1/4" surgical tubing  
Blue 438-1600 Red 438-1603  
Yellow 438-1601 Black 438-1604  
Green 438-1602



### Saliva Ejectors

(Unipack)

**Disposable Saliva Ejectors**  
Soft non-removable tip, vented tip  
100/pkg **\$2.95**  
White 436-1100  
Transparent 436-1101



(Winner)

1000/case **\$21.50**  
White/white 436-1002 Teal 436-1006  
Clear w/blue tips 436-1003 Raspberry 436-1007  
Green 436-1004 Sky Blue 436-1040  
Red 436-1005



(Starry Shine)

1000/case, 5 3/4" **\$21.50**  
White 436-1000  
Clear w/blue tips 436-1001

(Plasdent)

100/case, 5 3/4" **\$2.49**

#### Flexjet™

White, white tip 436-1010  
Clear, white tip 436-1011  
white, blue tip 436-1012  
Clear, blue tips 436-1013



#### Flexjet™ **\$2.49**

100/case, 5 3/4"  
Clear multicolor 436-1017  
Opaque multicolor 436-1018



# EVACUATION PRODUCTS

## Saliva Ejectors

(Plasdent)

### Flexform

100/case, 5 3/4"  
White w/clear adapter 436-0030 **\$4.95**



(Banta)

100/case, 5 3/4" **\$3.95**  
White 436-1015  
Clear w/blue tips 436-1016



(Crosstex)

Non-Removable Tip

**Comfort Plus™ Premium \$3.95**

100/bag  
White/White 436-0015  
Clear/Clear 436-0016  
Bubble Gum White/Pink 436-0017



**Comfort Plus™ Premium - Scented**

Bubble Gum White/Pink  
100/bag 436-0037 **\$3.95**



### Premium®

100/bag, 5 3/4" **\$3.95**  
Clear/Blue 436-0020  
Clear/White 436-0021  
White/White 436-0022  
Pink/Pink 436-0023  
Multicolor 436-0024



### Advantage

100/bag, 5 3/4" **\$3.49**  
Clear/Blue 436-0025  
Clear/Clear 436-0026  
White/White 436-0027  
Blue/Blue 436-0028



### Free Form / Hygoforming

Saliva ejector/tongue retractor  
White, 100/bag 436-0035 **\$7.95**  
Adapters, 10/bag 436-0036 **\$6.50**



### SE Cushion (Zirc)

Used in conjunction with the saliva ejector  
100/Pack 436-1200 **\$14.95**



### Lingua-Fix (Zirc)

Holds tongue in a steady position, helps  
protect the tongue from injuries  
10-Pack 436-1150 **\$9.95**  
48-Pack 436-1151 **\$39.50**



## Saliva Ejectors

### metal

(Vista)

**30° Angle \$21.95**

**Tapered Tip** (Reduces clogging)

**806001:** 1.5mm opening 436-1250  
Equivalent to 14P1-1/4

**806002:** 2.5mm opening 436-1251  
Equivalent to 14P2-1/4

**806003:** 3.0mm opening 436-1252  
Equivalent to 14P3-1/4



**90° Angle \$21.95**

**Tapered Tip** (Reduces clogging)

**806004:** 1.5mm opening 436-1253  
Equivalent to 14P1-90-1/4

**806005:** 2.5mm opening 436-1254  
Equivalent to 14P2-90-1/4

**806006:** 3.0mm opening 436-1255  
Equivalent to 14P3-90-1/4



## Surgical Tips

### plastic

(Winner)

Disposable, pkg of 25 **\$4.95**  
Small - White 1/8" 435-0001  
Large - Green 1/4" 435-0002  
large - Blue 1/4" 435-0003



(Unipack)

**Disposable Surgical Aspirator Tips**  
25/bag **\$4.95**

X-Small Blue 1/16" 435-0500  
Small White 1/8" 435-0501  
Large Green 1/4" 435-0502



## (StarryShine)

Disposable, pkg of 25	<b>\$4.95</b>
X-Small - Blue 1/16"	435-0005
Small - White 1/8"	435-0006
Large - Green 1/4"	435-0007
X-large - Blue 5/6"	435-0008



## (ibt Med)

Disposable, 25/pkg	<b>\$4.95</b>
Small White 1/8"	435-0008
Large Green 1/4"	435-0009

## (Pac-Dent)

Disposable or Autoclavable	<b>\$5.95</b>
25/pkg	
Green, large	435-0010
White, small	435-0011



## (Vista)

25/pkg	<b>\$6.95</b>
Small, White 1/8"	435-1010
Medium, Green 1/4"	435-1011
Large, Blue 3/8"	435-1012



## (Plasdent)

<b>Oralsurge II</b>	<b>\$6.95</b>
25/Bag	
<b>1/4" Orifice</b>	
Aquamarine	435-0021
White	435-0022
<b>1/8" Orifice</b>	
White	435-0023



## Oralsurge III

<b>3/8" Orifice</b>	<b>\$6.95</b>
vented tip, Seablue	435-0024

## Oralsurge-ST

<b>1/16" Orifice</b>	<b>\$8.95</b>
Blue	435-0025



## Oratip II (Pascal)

50/pkg	<b>\$12.95</b>
435-0026	

## Chromasurge™

<b>1/8" orifice, 100/Bag</b>	<b>\$11.60</b>
Assorted	435-0027




## Surgical Aspirators

metal

### (Vista)

Stainless steel, autoclavable  
Fits HVE Cut-Off Valve



### 3/16" Dia Tubing

**Tapered Tip (Reduces clogging) \$21.95**

**801001:** 1.5mm opening 437-0500  
Equivalent to 15P1A

**801002:** 2.5mm opening 437-0501  
Equivalent to 15P2A

**801003:** 3.0mm opening 437-0502  
Equivalent to 15P3A



**Slotted Tip (Reduces tissue-grab) \$21.95**

**801004:** 4.0mm opening 437-0505  
Equivalent to 16RSA

**Right Angled Tapered Tip \$21.95**

(Reduces clogging)  
**801005:** 1.5mm opening 437-0510  
Equivalent to 15P1-90A

**801006:** 2.5mm opening 437-0511  
Equivalent to 15P2-90A

**801007:** 3.0mm opening 437-0512  
Equivalent to 15P3-90A



**Right Angled Slotted Tip (Reduces tissue-grab)**

**801008:** 4.0mm opening 437-0515 **\$21.95**  
Equivalent to 16RS-90A

**1/4" Dia Tubing \$21.95**

**Tapered Tip (Reduces clogging)**

**801009:** 1.5mm opening 437-0520  
Equivalent to 46P1A

**801010:** 2.5mm opening 437-0521  
Equivalent to 46P2A

**801011:** 3.0mm opening 437-0522  
Equivalent to 46P3A



**5/6" Dia Tubing \$21.95**

**Tapered Tip (Reduces clogging)**

**801012:** 2.5mm opening 437-0530  
Equivalent to 56P2

**801013:** 3.0mm opening 437-0531  
Equivalent to 56P3

**801014:** 4.0mm opening 437-0532  
Equivalent to 56P4



**Frazier Style \$34.95**

**Fits 1/4" Surgical Tubing**

**805001:** 1.5mm opening 437-0535  
Equivalent to FR1

**805002:** 2.5mm opening 437-0536  
Equivalent to FR2

**805003:** 3.0mm opening 437-0537  
Equivalent to FR3



# EVACUATION PRODUCTS

## Evacuation System Cleaners

liquid

### E-Max (FCI)

Enzymatic super concentrated

1 Gallon 434-1030 **\$19.95**



### Pro E-Vac (Certol)

Concentrated

1 Gallon w/pump 434-1040 **\$34.99**

1oz Unit dose (24/box) 434-1041 **\$16.95**



### Vacusol Ultra (Biotrol)

Starter Kit 434-1044 **\$79.95**

Refill

32 oz Bottle 434-1045 **\$44.95**

96 oz Bottle 434-1046 **\$89.95**



### NeutraVac (Biotrol)

Starter Kit 434-1050 **\$69.50**

32 oz Bottle 434-1051 **\$39.95**

96 oz Bottle 434-1052 **\$99.95**



### PureVac (Sultan Chemists)

2 Liter Bottle 434-0055 **\$49.50**

5 Liter Bottle 434-0056 **\$99.95**

Walkabout Measuring & Dispensing Unit 434-0059 **\$69.95**



### PureVac Hg (Sultan Chemists)

1 oz makes 1 quart of solution

1 oz pump, non-Foaming & Biodegradable

21132 2-Liter (67oz) Bottle  
(67 treatments)

21135 5-Liter (169 oz) Bottle 434-0065 **\$49.95**

(169 treatments)  
434-0066 **\$99.95**



## Evacuation System Cleaners

powder

(California Dental)

2 lb powder 434-1000 **\$12.50**



### Citrizyme™ (Pascal)

Dual enzymatic cleaner for ultrasonic and evacuation systems, ultra concentrated

300g powder 434-1020 **\$27.50**

Unit Dose, 50/pkg 434-1025 **\$29.50**

900g powder 434-1026 **\$66.50**



## Evacuation Cups

### Liners for A-DEC

6oz white  
250/box 433-1000 **\$14.95**



### Dry Oral Cup (Plasdent) **\$2.95**

Beige, autoclavable, single pack

Type I - 8118 433-1010

Type II - 8119 433-1015

### Clear, autoclavable (Certol) **\$2.95**

Each 433-1020



## More Supplies For Less Money™





# EVACUATION PRODUCTS

## Disposable Traps

plastic

### ProEvac™ Disposable Traps (Pac-Dent)

144 filters/pack

**PD-5500** Dia 1 1/4" 438-1700  
fits A-Dec, Cascade & Performer

**PD-5501** Dia 2 1/8" 438-1701  
fits A-Dec w/plastic canisters

**PD-5502** Dia 2 1/4" 438-1702  
fits A-Dec w/metal canisters

**PD-5503** Dia 1 7/8" 438-1703  
fits Belmont, Biotec, Chayes,  
Den-Tal-Ez, Dentech, Pelton & Crane,  
Proma, etc

**PD-5505** Dia 3/4" 438-1704  
fits Pelton & Crane Spirit w/white  
plastic canisters

**PD-5508** Dia Universal 438-1705  
fits all saliva ejectors, 72/box

**PD-5509** fits DCI/Belmont 438-1706  
Realisty models, System 6000 models

**PD-5510** Dia 3" 438-1707  
fits Sirona - C2, C3, C4, C6, C8,  
Pelton & Crane S1, S2

**PD-5512** Dia 2 3/4" 438-1707  
fits A-dec 500

**PD-6100** Dia Universal 438-1707  
fits most Cuspidor Bowls  
including: Adec, Ampco, Belmont,  
Dabi Atlante, Dansereau, DCI,  
Den-Tal-Ez, Forest, Knight, Marus,  
MDT, Proma, Ritter, Schein

**PD-6200** Dia Universal 438-1707  
fits Siemens, Planmeca, and  
many cuspidors with small  
diameter openings

**PD-6500** Dia Universal 438-1707  
fits Kavo the environment

**\$19.95**



5500



5501



5502



5503



5505



5508



5509



5510



5512



6100



6200



6500

## Disposable Canisters

holder



### Disposable Canister (Plasdent) Each **\$8.95**

**DC8-2400**, fits Dentsply-MVS 0.75, 1.5, 4.0, 438-0030  
before 1991 (4 1/2" W x 5 1/8" H)

**DC8-2300**, fits Den-Tal-Ez dual after 1978, 438-0031  
Dentsply after 1991, Matrix 1000, 2000 &  
Dynamix models (3 1/2" W x 4 3/8" H)

**DC8-2350**, fits Apollo 2 hp, 3 hp single & 438-0032  
tandem, Air Techniques 1997 & later dual  
models, Vacudent-Workhorse dual & newer  
Vacstar (3 1/2" W x 4 3/8" H)

**DC8-2200**, fits Air Techniques-Vacstar, Apollo 438-0033  
1 hp 1.5 single, Tandem, 1, 2 hp Den-Tal-Ez  
single, Ohmeda-matrix, Turbine Industries,  
Techwest-Whirl Vac, Vacudent-Workhorse  
single, plus & other units. (2 3/4" W x 3 5/8" H)

### Easy-E-Trap (Pac-Dent) **\$5.95**

Kit consist of waste canister & lid, filter  
bowl, filter screen, replacement rubber  
gasket, & warning sticker

Model #2200 438-2000

Model #2300 438-2001

Model #2350 438-2002

Model #2400 438-2003

Model #2500 438-2004



### Mesh Filters (Pac-Dent)

Filters for central vacuum pump  
6 filters per pkg

**H: 3.75", dia 1.50"** 438-0050 **\$2.95**  
fits model #2200

**H: 4.25", dia 2.50"** 438-0051 **\$3.49**  
fits model #2300

**H: 4.25", dia 2.50"** 438-0052 **\$3.95**  
fits model #2350

