

Polishing Kits


Composite

Sof-lex Discs (3M ESPE)

1980: Contouring & Polishing Kit

120 discs in 1/2" and 3/8" diameters + one RA mandrel 440-1010 **\$99.50**

Refill Discs, 85/pkg **\$39.50**



1981	3/8"	440-1500	440-1501	440-1502	440-1503
1982	1/2"	440-1510	440-1511	440-1512	440-1513

Sof-lex XT Discs (3M ESPE)

2380: Contouring & Polishing Kit

120 discs in 1/2" and 3/8" diameters & one RA mandrel 440-1800 **\$99.50**

Refill Discs, 85/pkg **\$39.50**



2381	3/8"	440-1700	440-1701	440-1702	440-1703
2382	1/2"	440-1710	440-1711	440-1712	440-1713

Abrasive Paper Discs (Microdont)

Mixed Kit

40 Fine, 30 Medium & 30 Coarse 100/pkg

Small 10305.013	440-1600	\$12.95
Medium 10305.009	440-1630	\$13.95
Large 10305.008	440-1650	\$14.95



Refill Discs, 100/pkg

MONO FACE	Code#	Diameter
UNOFA	Code#	
	10305.014	15 mm
	10305.012	19 mm
	10305.002	22 mm
	10305.016	15 mm
	10305.011	19 mm
	10305.006	22 mm
	10305.015	15 mm
	10305.010	19 mm
	10305.004	22 mm

	Small	Medium	Large
	\$12.95	\$13.95	\$14.95
Fine	440-1601	440-1621	440-1641
	440-1602	440-1622	440-1642
	440-1603	440-1623	440-1643
Med	440-1604	440-1624	440-1644
	440-1605	440-1625	440-1645
	440-1606	440-1626	440-1646
Coarse	440-1607	440-1627	440-1647
	440-1608	440-1628	440-1648
	440-1609	440-1629	440-1649

Polident Kit (Microdont)

No metal centers, flexible Mylar discs, convenient Starter Kit, color coded strips, economically priced



Starter Kit 440-2000 **\$49.95**

84 Polishing discs, 35 Polishing strips, 1 silicone point, & 1 RA mandrel

Introductory Kit

28 pcs + mandrel 440-2010 **\$12.95**

Accessories, each

440-2020 **\$1.99**

10214.001

442-1105 **\$2.95**

10502.006



Accessories, Regular Discs

\$6.95

Refills, 24/box



	11001.001	● Top/X-Coarse	440-2050
	11001.002	● Bottom/X-Coarse	440-2051
	11001.003	● Top/Coarse	440-2052
	11001.004	● Bottom/Coarse	440-2053
	11001.005	● Top/Medium	440-2054
	11001.006	● Top/Fine	440-2055
	11001.007	● Top/X-Fine	440-2056

Polishing Discs (JSP)

Metal centers, flexible Mylar discs, convenient Starter Kit, color coded, 40 Polishing discs, 1 RA mandrel

440-2500 **\$19.95**



FINISHING & POLISHING

Polishing Kits composites


Diamond Polident Kit (Microdont)
Composite disc system. Thin mylar discs utilize diamond particles

Starter Kit 440-2300 **\$29.95**
21 pcs + mandrel



Snap-on Disc mandrel, 1/Pack

10502.006 442-1105 **\$2.99**



GRIT SIZE / GRANULACIÓN

- 25 μ Extra Fine Grain / Grano Extra Fino
- 40 μ Fine Grain / Grano Fino
- 90 ~ 125 μ Standard Grain / Grano Normal
- 125 ~ 150 μ Coarse Grain / Grano Grueso

Accessories, Regular Discs Refills,
24/box **\$13.95**

	11002.001	● Top/Coarse	440-2301
	11002.002	● Bottom/Coarse	440-2302
	11002.003	● Top/Medium	440-2303
	11002.004	● Bottom/Medium	440-2304
	11002.005	● Top/Fine	440-2305
	11002.006	● Bottom/Fine	440-2306
	11003.001	● Top/Coarse/Mini	440-2310
	11003.002	● Bottom/Coarse/Mini	440-2311
	11003.003	● Top/Medium	440-2312
	11003.004	● Bottom/Medium	440-2313
	11003.005	● Top/Fine	440-2314
	11003.006	● Bottom/Fine	440-2315

Poli Gloss (Microdont)

Snap-on silicon polishing system with three different grits

Starter Kit 440-2200 **\$17.95**
18 pcs + mandrel

Refills - 24/pkg **\$14.95**



KIT POLIGLOSS	Code#	Color
●	10217.003	440-2203
●	10217.004	440-2212
●	10217.006	440-2222
●	10217.007	440-2202
●	10217.008	440-2211
●	10217.010	440-2221

●	10217.011	440-2201
●	10217.012	440-2210
●	10217.014	440-2220

Speedy Gloss (Microdont)

Snap-on silicone point one-step system

Starter Kit 440-2100 **\$14.95**
18 pcs + mandrel

Refills **\$14.95**
24/pkg



KIT SPEEDYGLOSS	Code#	Color
●	10217.005	440-2110
●	10217.009	440-2111
●	10217.013	440-2112

Disc shape
Point shape
Cup shape



Venus® Supra (Heraeus)

Composite Polishing Kit
Venus Supra pre and high gloss polishers are made of silicone highly filled with microfine diamond powder. The five different shapes provide a universal solution for all anterior and posterior restorations.

Kit 440-2918 **\$79.50**
10 polishers in total: one of each shape in both pre-polishers & high gloss polishers.

Refill - RA Type **\$89.50**

Cup - Large	440-2920
Cup - Small	440-2921
Disc	440-2922
Flame - Large	440-2923
Flame - Small	440-2924



Polishing Strips

Polyester

Diamond Polyester (Microdont)

10/pkg **\$9.95**

	Narrow	Wide	Narrow	Wide	Narrow	Wide
# Microdont	10308.001	10308.002	10308.003	10308.004	10308.005	10308.006
Thickness Espesor (mm)	2,5	4,0	2,5	4,0	2,5	4,0
Length Comentarios/Comentarios	170	170	170	170	170	170
Grit Size Comentarios/Comentarios	20	20	40	40	60	60
Package Embalaje/Embalaje	10	10	10	10	10	10

Narrow

Green - Fine 440-3000
Blue - Medium 440-3002
Black - Coarse 440-3004

Wide

Green - Fine 440-3001
Blue - Medium 440-3003
Black - Coarse 440-3005

Abrasive Polyester (Microdont)

Pack of 100

Narrow **\$3.99**
Regular **\$4.99**
Wide **\$6.99**



MODEL Modelo	Code# Código	Thickness Espesor	Length Largo	Grit Size Comentarios	Packing Embalaje
	10302.001 10302.004	2.5 mm	170	240/320	Box/Caja 100 Box/Caja 150
	10302.002 10302.005	4.0 mm	170	240/320	Box/Caja 100 Box/Caja 150
	10302.003 10302.006	6.0 mm	170	240/320	Box/Caja 100 Box/Caja 150

Narrow 440-3020 Narrow 440-3030
Regular 440-3021 Regular 440-3031
Wide 440-3022 Wide 440-3032

Finishing Strips (Rite-Dent)

7" Aluminium Oxide

100/box

\$6.95

Narrow 440-4000

Medium 440-4010

Wide 440-4020



Sof-lex Finishing Strips (3M ESPE)

100/box

\$29.95

7" Aluminium Oxide

Fine/S-Fine #1956, 120/pkg 440-4100

Coarse/med #1954, 50/pkg 440-4101

Coarse/med #1954N, 100/pkg 440-4102



Spectra-A & Spectra-F Strips (Integra)

Spectra-A is made of resilient durable polyester fused with aluminum oxide grit, designed for use on the hardest composite materials. Spectra-F is used for general purpose finishing & polishing. Effective wet or dry, full 1/2" open center, color coded

Spectra-A, box of 150

\$17.95

Fine 440-2650

Medium 440-2651

Coarse 440-2652

Assorted - Medium width 440-2653

Spectra-F, box of 150

\$17.95

Fine 440-2660

Medium 440-2661

Coarse 440-2662

Assorted - Medium width 440-2663



Dual Alo #54

\$21.50

Evenwet® Strips (Integra)

7" long w medium grit on one half and fine grit on the other half. Effective wet or dry. 150 strips / box

Narrow, 2.5mm 440-2700

Medium 3.18mm 440-2701

Wide 4.00mm 440-2702



Diamond Blade w/Saw Blade

Diamond Strips (Edenta)

for approximal contouring of composite, acrylic, glass ionomer cement & amalgams
3.75mm width with 40mm of sawblade

10/pkg 440-2710 **\$41.95**



FINISHING & POLISHING

Polishing Strips

metal

Stainless Steel Diamond Strip (Microdont)

10/pkg **\$19.95**



MODEL	Code#	Thickness	Length	Grit Size	Packing
MODEL	Code#	mm	Inches	Conclusion	Strips
10309.003 10309.005 10309.007		2.5 mm	150	(20) Fine/Fino (40) Medium/Medio (80) Coarse/Gruoso	Box/Caja 10
10309.004 10309.006 10309.008		4.0 mm	150	(20) Fine/Fino (40) Medium/Medio (80) Coarse/Gruoso	Box/Caja 10

Narrow

Wide

Fine 440-3060

440-3064

Medium 440-3061

440-3065

Coarse 440-3062

440-3066

Stainless Steel Super Diamond Strip (Microdont)

10/pkg **\$19.95**

SUPER DIAMOND STRIPS - SUPER TIRAS DE ACERO DIAMANTADA					
MODEL	Code#	Thickness	Length	Grit Size	Packing
MODEL	Code#	mm	Inches	Conclusion	Strips
10309.001		2.5 mm	150	45	Box/Caja 10
10309.002		4.0 mm	150	45	Box/Caja 10

Narrow

Wide

Medium 440-3063

440-3067

Stainless Steel Abrasive Strip (Microdont)

Extra Coarse, Pack of 12



MODEL	Code#	Thickness	Length	Grit Size	Packing
MODEL	Code#	mm	Inches	Conclusion	Strips
MONO FACE - MONO CARA					
10304.004		2.5 mm	150	240	Box/Caja 12
10304.005		4.0 mm	150	240	Box/Caja 12
10304.006		6.0 mm	150	240	Box/Caja 12
DOUBLE FACE - DOBLE CARA					
10304.010		2.5 mm	150	240	Box/Caja 12
10304.008		4.0 mm	150	240	Box/Caja 12
10304.009		6.0 mm	150	240	Box/Caja 12

Monoface \$5.95

Double face \$8.95

Narrow 440-3070

Narrow 440-3075

Regular 440-3071

Regular 440-3076

Wide 440-3072

Wide 440-3077

Metal Strips (IDS)

One sided, 12/box **\$9.95**

4mm 440-7200

6mm 440-7210



Abrasive Silicone

amalgam

Polishing Kit - RA Type (Microdont)

Clear plastic box

Complete Kit

14/pack 440-2410 **\$14.95**

Starter Kit

12/pack 440-2412 **\$14.95**

Intro Kit

8/pack 440-2414 **\$12.95**

Accessories Refills

12/pack **\$14.95**

Grit Size / Granulación / Granulometria

- 1200 Mesh / Extra Fine Grain / Grano Extra Fino / Grão Extra Fino
- 1000 Mesh / Fine Grain / Grano Fino / Grão Fino
- 800 Mesh / Standard Grain / Grano Normal / Grão Normal
- 400 Mesh / Coarse Grain / Grano Gruoso / Grão Grosso



POINTS	Code#	Grit Size
POINTS	Code#	Granulación

	10204.003	●
	10204.004	●
	10204.001	●
	10204.002	●

Polishing Points

RA	FG
442-1000	442-1004
442-1001	442-1005
442-1002	442-1006
442-1003	442-1007

CUP	Code#	Grit Size
CUP	Code#	Granulación

	10203.003	●
	10203.004	●
	10203.001	●
	10203.002	●

Polishing Cups

RA	FG
442-1010	442-1014
442-1011	442-1015
442-1012	442-1016
442-1013	442-1017

KNIFE EDGE	Code#	Grit Size
KNIFE	Code#	Granulación

	10202.003	●
	10202.004	●
	10202.001	●
	10202.002	●

Polishing Discs

RA	FG
442-1030	442-1034
442-1031	442-1035
442-1032	442-1036
442-1033	442-1037

FINISHING & POLISHING

Abrasive Silicone

amalgam

Brownie®, Greenie®, Supergreenie (Shofu)

12/pkg

\$19.65

	Cup	Mini Point	KN7	PC2
Brownie® - Pre-Polish				
RA	440-5000	440-5001	440-5002	440-5003
FG	440-5005	440-5006	-	-
Greenie® - Polish				
RA	440-5020	440-5021	440-5022	440-5023
FG	440-5025	440-5026	-	-
Supergreenie® - Super Polish				
RA	440-5050	440-5051		
FG	440-5055	440-5056		

KNIFE EDGE LEPWA	Code# ClAppt	Grit Size Granulatin
	10202.009	●
	10202.010	●
	10202.011	●

Polishing Discs 12/pack

\$14.95

RA	FG
440-2480	440-2483
440-2481	440-2484
440-2482	440-2485

Polishing Stone

white

Starter Kit (Microdont)

Clear plastic box, 12/pack

440-2450

\$14.95

Aluminum box, 12/pack

440-2455

\$19.95



Coarse Medium Fine



POINTS PUNTS	Code# ClAppt	Grit Size Granulatin
	10204.009	●
	10204.010	●
	10204.011	●

Polishing Points Refill, 12/pack

\$14.95

RA	FG
440-2460	440-2463
440-2461	440-2464
440-2462	440-2465

CUP COPA	Code# ClAppt	Grit Size Granulatin
	10203.009	●
	10203.010	●
	10203.011	●

Polishing Cups Refill, 12/pack

\$14.95

440-2470	440-2473
440-2471	440-2474
440-2472	440-2475

Dura-White (Shofu)

12/pkg

\$17.95

Shape		CN1		CY2		FL1
FG		440-5200		440-5202		-
RA		440-5201		-		440-5203
Shape		FL2		FL3		IC2
FG		440-5220		-		440-5223
RA		440-5221		440-5222		-
Shape		IC4		RD1		RD2
FG		-		440-5231		440-5233
RA		440-5230		440-5232		440-5234
Shape		RE1		TC1		TC4
FG		-		440-5251		-
RA		440-5250		440-5252		440-5253

Polishing Stone

green

Dura-Green Stone (Shofu)

12/pkg

\$17.95

Shape		CN1		CY2		FL1
FG		440-5400		440-5402		-
RA		440-5401		-		440-5404
Shape		FL2		FL3		RD1
FG		440-5420		-		440-5423
RA		440-5421		440-5422		440-5424
Shape		RD2		RE1		TC1
FG		440-5425		-		440-5428
RA		440-5426		440-5427		440-5429



FINISHING & POLISHING

Abrasive Silicone

porcelain

CeraMaster® (Shofu)

Diamond impregnated polishers, for polishing porcelain crowns, 3/pack



Kits

Assorted Fine **RA**
\$39.95
Bullet, Knife, Cup, Mini Point 440-2815

Assorted Coarse
Bullet, Knife, Cup, Mini Point 440-2816

Refills, each shape **\$29.95**
3/pack




	Bullet	Knife	Cup	Mini-Pt
ISO#	040	065	090	030
Fine				
CA/RA	440-2801	440-2802	440-2803	440-2804
Coarse				
CA/RA	440-2821	440-2822	440-2823	440-2824

Abrasives

kits

CADMaster CAD/CAM™ (Shofu)

FG Kit 440-2900 **\$89.50**
Dura-green stones (FL2 & 2 x RD2),
Dura-white stones (FL2 & 2 x RD2),
Ceramaster (2 x mini points),
Ceramaster Coarse, Diamond Robot
Point & T&F Hybrid Point (7406)



RA Kit 440-2902 **\$109.95**
Dura-green stones (FL2, RD2, TC4 & WH6),
Dura-white stones (FL2, RD2 & TC4),
Ceramaster (bullet & WH6), Ceramaster
Coarse (bullet & WH6)

HP Kit 440-2904 **\$79.95**
Dura-green stones (FL2 & 2 x RD2), Dura-white
stones (FL2 & 2 x RD2), Ceramaster (bullet,
knife, cup & mini-point), Ceramaster Coarse
(bullet, knife, cup & mini-point)



CompoMaster® Finishing & Polishing Kit (Shofu)

FG hybrid T&F points (7106, 7406, 7901, LT1),
compomaster coarse polishing points (knife,
cup & mini-point), compomaster polishing
points (knife, cup & mini-point)



0331 Classic Plastic Kit - CA/FG Shank 440-2914 **\$39.95**

0331B Aluminum bur Kit - CA/FG 440-2915 **\$39.95**

RA Mandrels (Microdont)

RA	Code/Color	ISO	Ø	Length	Description
	10502.001	611 ~ [323]	023	13,0	FG Converted Adapter FG
	10502.002	603 ~ [303]	050	22,0	Discs and Saws Discs & Saws
	10502.003	611 ~ [322]	050	31,0	Polisher and Felts Polisher & Felts
	10502.004	603 ~ [303]	030	19,0	Discs and Saws - Mini Discs & Saws - Mini
	10502.005	615 ~ [313]	035	25,0	Discs - Pop-On Discs - Pop-On
	10502.006	611 ~ [327]	045	23,0	Special Discs - Snap Discs Expandable - Snap
	10502.007	611 ~ [327]	045	45,0	Special Discs - Snap Long Discs Expandable - Snap Long
	10502.009	612 ~ [327]	050	23,0	Occlusal Polisher Fusion Chisels

Disc Saw, 1/Pack 442-1100 **\$0.99**
Disc Saw, 1/Pack 442-1101 **\$0.99**
Pop-on Disc, 3/Pack 442-1102 **\$14.99**
Occlusal Polisher, 1/Pack 442-1103 **\$0.99**
Polisher Felt, 1/Pack 442-1104 **\$0.99**
Snap-on Disc, 1/Pack 442-1105 **\$0.99**
Snap-on Disc, 3/Pack 442-1106 **\$5.99**
FG Converted, 1/Pack 442-1107 **\$4.99**

Composite Polishing

interproximal

CPR Stain Remover (ICCare)

Perfect for composites, veneers
and restorations, non-pumice
100g cups/jar

Mint 440-5800 **\$59.95**



DiaGlaze (Deldent)

Assorted Diamond polishing paste, 3 x 2g
syringe. Extra Fine, Fine, Regular

Kit 440-5810 **\$47.95**

Each Syringe, 2g **\$19.95**

X-Fine 440-5811

Fine 440-5812

Regular 440-5813



Disc & Separator


interproximal

Interproximal Contact Discs

(Link Industries)

Trial Kit \$21.95
1 of each Smooth, Smooth/Abrasive, and Abrasive 442-1275

Refill Discs, 8/pkg
Smooth 442-1276 **\$33.50**
Smooth/Abrasive 442-1277 **\$57.50**
Abrasive 442-1278 **\$49.95**



Proxidisc (Centrix) \$66.95

Proximal separator discs, 18/pkg

Smooth 440-2750
Abrasive 440-2751



Cerisaw (Denmat)

Proximal separator saw, single

Anterior kit with blade 440-2755 **\$39.95**
#031336510

Posterior kit with blade 440-2756 **\$39.95**
#031336610

Anterior blade, 10 pc 440-2757 **\$36.95**
#031336520

Posterior blade, 10 pc 440-2758 **\$36.95**
#13371



For the latest in Composite Finishing & Polishing

It's Your Choice.



Speedy

Economical, fast and easy snap-on silicone point system that produces a satin-like polish with one point

AVAILABLE IN 3 SNAP-ON SHAPES
CONVENIENT STARTER KIT AVAILABLE
REFILLS AVAILABLE IN PACKS OF 24

18 PIECES
+
MANDREL



Poli

Utilizes a 3-grit sequence that works great on any composite. The combination produces superior esthetic results!

REFILLS AVAILABLE IN PACKS OF 24

18 PIECES
+
MANDREL



the Latest in Snap-On Silicone Systems



POLIDONT

DIAMOND DISC SYSTEM

The latest in Composite Disk Systems

UTILIZES DIAMOND PARTICLES TO ACHIEVE THE ULTIMATE RESULTS

THIN MYLAR DISCS PLUS MINI DISCS TO REACH DIFFICULT AREAS

3 EASY STEPS ONLY 3 COLOR-CODED GRITS TO ACHIEVE BEST RESULTS

CONVENIENT STARTER KIT AVAILABLE

REFILLS AVAILABLE IN PACKS OF 24



MICRODONT USA



28 PIECES
+
MANDREL