

## Polishing Kits

### Composite

#### Sof-lex Discs (3M ESPE)

##### #1980: Contouring & Polishing Kit

120 discs in 1/2" and 3/8" diameters + one RA mandrel 440-1010 \$162.95

Refill Discs, 85/pkg \$65.95



1981	3/8"	440-1500	440-1501	440-1502	440-1503
1982	1/2"	440-1510	440-1511	440-1512	440-1513

#### Sof-lex XT Discs (3M ESPE)

##### #2380: Contouring & Polishing Kit

120 discs in 1/2" and 3/8" diameters & one RA mandrel 440-1800 \$162.95

Refill Discs, 85/pkg \$65.95



2381	3/8"	440-1700	440-1701	440-1702	440-1703
2382	1/2"	440-1710	440-1711	440-1712	440-1713

#### Sof-lex Spiral (3M ESPE)

##### Spiral Finishing & Polishing Wheels

##### #5082SK: System Kit, 10 Sof-Lex Spirala

Finishing Wheels, 10 Sof-Lex Spiral Polishing Wheels, 30 Each Sof-Lex XT Polishing Discs 12.7mm and 9.6mm (Coarse, Medium, Fine, Superfine), 15 Sof-Lex Finishing Strips, 1 RA Mandrel



440-1900 \$199.95

#5082I: Intro Kit, 10 Sof-Lex 440-1901 \$53.95

Spiral Finishing Wheels, 10 Sof-Lex Spiral Polishing Wheels, 1 RA Mandrel, 1 IFU Spiral Wheels, 1 IFU Sof-Lex Mandrel

#5080: Finishing Refill, 30/Pkg 440-1910 \$66.95

#5081: Polishing Refill, 30/Pkg 440-1911 \$66.95

#### SUPER-SNAP SINGLES (Shofu)

Refill Kit 48/pkg \$159.95

Standard PN 0502 440-2070

Singles Mini PN 0502 440-2071



#### SUPER-SNAP Disk Replacements (Shofu)

Refill Kit, 50/pkg \$23.95

- Black/Violet: Silicon Carbide grit
- Green/Red: Aluminum Oxide grit

SHAPE CATEGORY	Contouring		Finishing				Polishing		Super-Polishing		
SHAPE											
COLOR/GRIT	Black (coarse)		Violet (medium)				Green (fine)		Red (super fine)		
REFILLS (50 PK)	PN	L506	L507	L508	L528	L509	L522	L501	L519	L502	L521

#L506	440-2080
#L507	440-2081
#L508	440-2082
#L528	440-2083
#L509	440-2084
#L522	440-2085
#L501	440-2086
#L519	440-2087
#L502	440-2088
#L521	440-2089

#### SUPER-SNAP Polystrips (Shofu)

Refill Kit, 100/pkg \$26.95

- Black/Violet: coarse/medium
- Green/Red: fine/super fine

SHAPE CATEGORY	Contouring	Finishing	Polishing	Super-Polishing
SHAPE				
COLOR/GRIT	Black/Violet (coarse/medium)		Green/Red (fine/super fine)	
(100 PK)	PN	L525		L526

#### SUPER-SNAP SuperBuff Disks (Shofu)

\$27.95

#0535	25 Disks, 2 CA Mandrels	440-2030
#0536	50 Disks	440-2031
#0439	6 CA Metal Mandrel	440-2032
#0440	100 CA Plastic Mandrel	440-2033



# FINISHING & POLISHING

## Polishing Kits

### Composite

**SUPER-SNAP Buff Disks (Shofu)**  
For composites and porcelain polishing



#0523	25 Disks, Regular	440-2040	\$36.95
#0536	25 Disks, Mini	440-2041	\$41.95

**OneGloss (Shofu)**  
Composites, temporary crowns, removal of orthodontic cements and removal of stains from enamel



#0180	25 Cup, 25 Mini Point, 25 IC shapes, 3 CA Mandrels	440-2045	\$84.95
-------	--	----------	---------

**OneGloss PS (Shofu)**  
Ideal finish to all composite restorations with a plastic shank.



#0175	15 Cup, 20 Mini Point, 15 IC shapes	440-2046	\$104.95
-------	-------------------------------------	----------	----------

SHAPE CATEGORY	Cup	Midi-Pt.	IC	Mini-Pt.	Assort. Kit	CA Mandrel
SHAPE					See Descriptions Above	
ISO #	060	060	100	030		
<b>OneGloss</b>						
CA (50 pk)	PN 0181	0182	0183	0186 (12 pk)	0180	0435 (6 pk)
<b>OneGloss PS</b>						
CA (50 pk)	PN 0176	0177	0178	-	0175	-

**OneGloss - Refill, 50/Pkg (Shofu)** \$69.95

#0181	440-2090
#0182	440-2091
#0183	440-2092
#0186	440-2093
#0435	440-2094

**OneGloss PS- Refill, 50/Pkg (Shofu)** \$104.95

#0176	440-2060
#0177	440-2061
#0178	440-2062

## Abrasive Paper Discs (Microdont)

Mixed Kit  
40 Fine, 30 Medium & 30 Coarse  
100/pkg

Small 10305.013	440-1600	\$21.95
Medium 10305.009	440-1630	\$21.95
Large 10305.008	440-1650	\$23.95



## Refill Discs, 100/pkg

MONO FACE	Color	Diameter
UNI-DENT	Blue	
	10305.014	15 mm
	10305.012	19 mm
	10305.002	22 mm
	10305.016	15 mm
	10305.011	19 mm
	10305.006	22 mm
	10305.015	15 mm
	10305.010	19 mm
	10305.004	22 mm

	Small	Medium	Large
Fine	440-1601	440-1621	440-1641
	440-1602	440-1622	440-1642
	440-1603	440-1623	440-1643
Med	440-1604	440-1624	440-1644
	440-1605	440-1625	440-1645
	440-1606	440-1626	440-1646
Coarse	440-1607	440-1627	440-1647
	440-1608	440-1628	440-1648
	440-1609	440-1629	440-1649

**Introductory Kit**  
28 pcs + mandrel 440-2010

Accessories, each	
	440-2020
10214.001	
	440-2025
10502.006	



## Accessories, Regular Discs

Refills, 24/box

	11001.001	● Top/X-Coarse	440-2050
	11001.002	● Bottom/X-Coarse	440-2051
	11001.003	● Top/Coarse	440-2052
	11001.004	● Bottom/Coarse	440-2053
	11001.005	● Top/Medium	440-2054
	11001.006	● Top/Fine	440-2055
	11001.007	● Top/X-Fine	440-2056



# FINISHING & POLISHING

**Diamond Polident Kit (Microdont)**  
Composite disc system. Thin mylar discs utilize diamond particles

Starter Kit 440-2300 \$35.95  
21 pcs + mandrel



Snap-on Disc mandrel, 1/Pack

10502.006 442-1105 \$5.95



### GRIT SIZE / GRANULACIÓN

- 25 µ Extra Fine Grain / Grano Extra Fino
- 40 µ Fine Grain / Grano Fino
- 90 - 125 µ Standard Grain / Grano Normal
- 125 - 150 µ Coarse Grain / Grano Grueso

Accessories, Regular Discs Refills,  
24/box \$19.95

	11002.001	● Top/Coarse	440-2301
	11002.002	● Bottom/Coarse	440-2302
	11002.003	● Top/Medium	440-2303
	11002.004	● Bottom/Medium	440-2304
	11002.005	● Top/Fine	440-2305
	11002.006	● Bottom/Fine	440-2306
	11003.001	● Top/Coarse/Mini	440-2310
	11003.002	● Bottom/Coarse/Mini	440-2311
	11003.003	● Top/Medium	440-2312
	11003.004	● Bottom/Medium	440-2313
	11003.005	● Top/Fine	440-2314
	11003.006	● Bottom/Fine	440-2315

**Polident Kit (Microdont)**  
No metal centers, flexible Mylar discs, convenient Starter Kit, color coded strips, economically priced

Starter Kit 440-2000 \$69.95  
84 Polishing discs, 35 Polishing strips,  
1 silicone point, & 1 RA mandrel



**Poli Gloss (Microdont)**  
Snap-on silicon polishing system with three different grits

Starter Kit 440-2200 \$19.95  
18 pcs + mandrel

Refills - 24/pkg \$19.95

REFILLS / PZS	Color	Code
●		10217.003
●		10217.004
●		10217.005
●		10217.007
●		10217.008
●		10217.010

440-2203

440-2212

440-2222

440-2202

440-2211

440-2221

●		10217.011	440-2201
●		10217.012	440-2210
●		10217.014	440-2220



**Speedy Gloss (Microdont)**

Snap-on silicone point one-step system

Starter Kit 440-2100 \$19.95  
18 pcs + mandrel

Refills 24/pkg \$19.95

REFILLS / PZS	Color	Code
●		10217.006
●		10217.009
●		10217.013

440-2110

440-2111

440-2112

Disc shape  
Point shape  
Cup shape



**Venus® Supra (Heraeus)**  
Composite Polishing Kit

Venus Supra pre and high gloss polishers are made of silicone highly filled with microfine diamond powder. The five different shapes provide a universal solution for all anterior and posterior restorations.

Kit 440-2918 \$109.95

10 polishers in total: one of each shape in both pre-polishers & high gloss polishers.

Refill - RA Type \$114.95

Cup - Large 440-2920

Cup - Small 440-2921

Disc 440-2922

Flame - Large 440-2923

Flame - Small 440-2924




# FINISHING & POLISHING

## Polishing Strips

Polyester

### Diamond Polyester (Microdont)

10/pkg \$12.95



	Narrow	Wide	Narrow	Wide	Narrow	Wide
# Microdont	10308.001	10308.002	10308.003	10308.004	10308.005	10308.006
Thickness (Espesor / Espaisseur)	2,5	4,0	2,5	4,0	2,5	4,0
Length (Comprimento / Longueur)	170	170	170	170	170	170
Grit Size (Granulometría / Granulométrie)	20	20	40	40	60	60
Package (Embalaje / Emballage)	10	10	10	10	10	10

Narrow		Wide	
Green - Fine	440-3000	Green - Fine	440-3001
Blue - Medium	440-3002	Blue - Medium	440-3003
Black - Coarse	440-3004	Black - Coarse	440-3005

### Abrasive Polyester (Microdont)

Pack of 100

Narrow	\$11.95
Regular	\$11.95
Wide	\$12.95



MODEL	Color / Color	Thickness / Espesor	Length / Largo	Grit Size / Granulometría	Package / Paquete
10302.001 10302.004	Green	2.5 mm	170	240/320	Box/Caja 100 Box/Caja 150
10302.002 10302.005	Blue	4.0 mm	170	240/320	Box/Caja 100 Box/Caja 150
10302.003 10302.006	Black	8.0 mm	170	240/320	Box/Caja 100 Box/Caja 150

Narrow	440-3020	Narrow	440-3030
Regular	440-3021	Regular	440-3031
Wide	440-3022	Wide	440-3032

### Finishing Strips (Rite-Dent)

7" Aluminum Oxide

100/box \$6.95

Narrow	440-4000
Medium	440-4010
Wide	440-4020



### Sof-Iex Finishing Strips (3M ESPE)

100/box \$43.95

7" Aluminium Oxide

Fine/S-Fine #1956, 120/pkg	440-4100
Coarse/med #1954, 50/pkg	440-4101
Coarse/med #1954N, 100/pkg	440-4102



### Spectra-A & Spectra-F Strips (Integra)

Spectra-A is made of resilient durable polyester fused with aluminum oxide grit, designed for use on the hardest composite materials. Spectra-F is used for general purpose finishing & polishing. Effective wet or dry, full 1/2" open center, color coded

Spectra-A, box of 150 \$22.95

Fine	440-2650
Medium	440-2651
Coarse	440-2652
Assorted - Medium width	440-2653

Spectra-F, box of 150 \$22.95

Fine	440-2660
Medium	440-2661
Coarse	440-2662
Assorted - Medium width	440-2663



Dual Alo #54 \$22.95

Evenwet® Strips (Integra)

7" long w medium grit on one half and fine grit on the other half. Effective wet or dry. 150 strips / box

Narrow, 2.5mm	440-2700
Medium 3.18mm	440-2701
Wide 4.00mm	440-2702



## Diamond Blade w/Saw Blade

### Diamond Strips (Edenta)

for approximal contouring of composite, acrylic, glass ionomer cement & amalgams  
3.75mm width with 40mm of sawblade

10/pkg 440-2710 \$50.95





# FINISHING & POLISHING

## Polishing Strips

metal

### Stainless Steel Diamond Strip (Microdont)

10/pkg \$19.95



MODEL	Code#	Thickness	Length	Grit Size	Packing
	10309.003	2.5 mm	150	200 Fine/Fino	Box/Caja 10
	10309.005				
	10309.007				
	10309.004	4.0 mm	150	140 Fine/Finco	Box/Caja 10
	10309.006				
	10309.008				

Narrow

Wide

Fine	440-3060	440-3064
Medium	440-3061	440-3065
Coarse	440-3062	440-3066

### Stainless Steel Super Diamond Strip (Microdont)

10/pkg \$19.95

MODEL	Code#	Thickness	Length	Grit Size	Packing
	10309.001	2.5 mm	150	45	Box/Caja 10
	10309.002	4.0 mm	150	45	Box/Caja 10

Narrow

Wide

Medium	440-3063	440-3067
--------	----------	----------

### Stainless Steel Abrasive Strip (Microdont)

Extra Coarse, Pack of 12



MODEL	Code#	Thickness	Length	Grit Size	Packing
<b>MONO FACE - MONO CARA</b>					
	10304.004	2.5 mm	150	240	Box/Caja 12
	10304.005	4.0 mm	150	240	Box/Caja 12
	10304.006	6.0 mm	150	240	Box/Caja 12
<b>DOUBLE FACE - DOBLE CARA</b>					
	10304.010	2.5 mm	150	240	Box/Caja 12
	10304.008	4.0 mm	150	240	Box/Caja 12
	10304.009	6.0 mm	150	240	Box/Caja 12

Monoface	\$5.95	Double face	\$8.95
Narrow	440-3070	Narrow	440-3075
Regular	440-3071	Regular	440-3076
Wide	440-3072	Wide	440-3077

### Metal Strips (IDS)

One sided, 12/box \$9.95

4mm	440-7200
6mm	440-7210



## Abrasive Silicone

amalgam

### Polishing Kit - RA Type (Microdont)

Clear plastic box

#### Complete Kit

14/pack 440-2410 \$14.95

#### Starter Kit

12/pack 440-2412 \$14.95

#### Intro Kit

8/pack 440-2414 \$14.95

#### Accessories Refills

12/pack

\$15.95

Grit Size / Granulación / Granulometría

- 1200 Mesh / Extra Fine Grain / Grano Extra Fino / Grão Extra Fino
- 1000 Mesh / Fine Grain / Grano Fino / Grão Fino
- 800 Mesh / Standard Grain / Grano Normal / Grão Normal
- 400 Mesh / Coarse Grain / Grano Grosso / Grão Grosso

POINTS	Code#	Grit Size
--------	-------	-----------

POINTS	Code#	Grit Size
--------	-------	-----------

POINTS	Code#	Grit Size
--------	-------	-----------

POINTS	Code#	Grit Size
--------	-------	-----------

POINTS	Code#	Grit Size
--------	-------	-----------

CUP	Code#	Grit Size
-----	-------	-----------

CUP	Code#	Grit Size
-----	-------	-----------

CUP	Code#	Grit Size
-----	-------	-----------

CUP	Code#	Grit Size
-----	-------	-----------

CUP	Code#	Grit Size
-----	-------	-----------

KNIFE EDGE	Code#	Grit Size
------------	-------	-----------

KNIFE EDGE	Code#	Grit Size
------------	-------	-----------

KNIFE EDGE	Code#	Grit Size
------------	-------	-----------

KNIFE EDGE	Code#	Grit Size
------------	-------	-----------

#### Polishing Points

RA	FG
----	----

442-1000 442-1004

442-1001 442-1005

442-1002 442-1006

442-1003 442-1007

442-1003 442-1007

#### Polishing Cups

RA	FG
----	----

442-1010 442-1014

442-1011 442-1015

442-1012 442-1016

442-1013 442-1017

442-1013 442-1017

#### Polishing Discs

RA	FG
----	----

442-1030 442-1034

442-1031 442-1035

442-1032 442-1036

442-1032 442-1036

442-1033 442-1037



# FINISHING & POLISHING

## Polishing Stone

FG

### Allbeing International

#### Dura-White

12/pkg \$14.95

Shape	CN1	CY2	FL1
FG	440-5300	440-5302	-
RA	440-5301	-	440-5303
Shape	FL2	FL3	IC2
FG	440-5320	-	440-5323
RA	440-5321	440-5322	-
Shape	IC4	RD1	RD2
FG	-	440-5331	440-5333
RA	440-5330	440-5332	440-5334
Shape	RE1	TC1	TC4
FG	-	440-5351	-
RA	440-5350	440-5352	440-5353

#### Dura-Green Stone (Allbeing)

12/pkg \$14.95

Shape	CN1	CY2	FL1
FG	440-5500	440-5502	-
RA	440-5501	-	440-5504
Shape	FL2	FL3	RD1
FG	440-5520	-	440-5523
RA	440-5521	440-5522	440-5524
Shape	RD2	RE1	TC1
FG	440-5525	-	440-5528
RA	440-5526	440-5527	440-5529

#### Starter Kit (Microdont)

Clear plastic box, 12/pack  
440-2450

Aluminum box, 12/pack  
440-2455



Coarse Medium Fine  
● ● ●

POINTS	Code#	Grit Size
POINT	Color	Grain Size
	10204.009	●
	10204.010	●
	10204.011	●

#### Polishing Points Refill, 12/pack

RA	FG
440-2460	440-2463
440-2461	440-2464
440-2462	440-2465

CUP	Code#	Grit Size
CUP	Color	Grain Size
	10203.009	●
	10203.010	●
	10203.011	●

#### Polishing Cups Refill, 12/pack

440-2470	440-2473
440-2471	440-2474
440-2472	440-2475

EDGE	Code#	Grit Size
EDGE	Color	Grain Size
	10202.009	●
	10202.010	●
	10202.011	●

#### Polishing Discs 12/pack

RA	FG
440-2480	440-2483
440-2481	440-2484
440-2482	440-2485



# FINISHING & POLISHING

## Abrasive Silicone

amalgam

Brownie®, Greenie®, Supergreenie (Shofu)

12/pkg

\$29.95

	Cup	Mini Point	KN7	PC2
<b>Brownie® - Pre-Polish</b>				
RA	440-5000	440-5001	440-5002	440-5003
FG	440-5005	440-5006	-	-
<b>Greenie® - Polish</b>				
RA	440-5020	440-5021	440-5022	440-5023
FG	440-5025	440-5026	-	-
<b>Supergreenie® - Super Polish</b>				
RA	440-5050	440-5051		
FG	440-5055	440-5056		

## Polishing Stone

green

Dura-Green Stone (Shofu)

12/pkg

\$25.95

Shape	CN1	CY2	FL1
FG	440-5400	440-5402	-
RA	440-5401	-	440-5404
Shape	FL2	FL3	RD1
FG	440-5420	-	440-5423
RA	440-5421	440-5422	440-5424
Shape	RD2	RE1	TC1
FG	440-5425	-	440-5428
RA	440-5426	440-5427	440-5429

## Polishing Stone

white

Dura-White (Shofu)

12/pkg

\$25.95

Shape	CN1	CY2	FL1
FG	440-5200	440-5202	-
RA	440-5201	-	440-5203
Shape	FL2	FL3	IC2
FG	440-5220	-	440-5223
RA	440-5221	440-5222	-
Shape	IC4	RD1	RD2
FG	-	440-5231	440-5233
RA	440-5230	440-5232	440-5234
Shape	RE1	TC1	TC4
FG	-	440-5251	-
RA	440-5250	440-5252	440-5253

## Disc & Separator

interproximal

Interproximal Contact Discs

(Link Industries)

Trial Kit

\$23.95

1 of each Smooth, Smooth/Abrasive, and Abrasive


442-1275

Refill Discs, 8/pkg

Abrasive

442-1278

\$66.95



COUTURE™ PROXIDISCS® (Centrix)

Proximal separator discs, 18/pkg

Smooth 2-sided, 8/pkg 440-2750 \$56.95

Diamond 1-sided, 4/pkg 440-2751 \$82.95



Cerisaw (Denmat)

Proximal separator saw, single

Anterior kit with blade 440-2755 \$49.95

#031336510

Posterior kit with blade 440-2756 \$49.95

#031336610

Anterior blade, 10 pc 440-2757 \$45.95

#031336520

Posterior blade, 10 pc 440-2758 \$45.95

#13371



# FINISHING & POLISHING

## Abrasive Silicone

porcelain






### CeraMaster® (Shofu)

Diamond impregnated polishers, for polishing porcelain crowns, 3/pack




Kits

	RA
Assorted Fine Bullet, Knife, Cup, Mini Point	\$58.95 440-2815
Assorted Coarse Bullet, Knife, Cup, Mini Point	440-2816
Refills, each shape 3/pack	\$43.95

	Bullet	Knife	Cup	Mini-Pt
fine				
coarse				
ISO#	040	065	090	030
<b>Fine</b>				
CA/RA	440-2801	440-2802	440-2803	440-2804
<b>Coarse</b>				
CA/RA	440-2821	440-2822	440-2823	440-2824

### CompoMaster® Finishing & Polishing Kit (Shofu)

FG hybrid T&F points (7106, 7406, 7901, LT1), compomaster coarse polishing points (knife, cup & mini-point), compomaster polishing points (knife, cup & mini-point) - 10 burs



0331 Classic Plastic Kit - CA/FG Shank	440-2914	\$39.95
0331B Aluminum bur Kit - CA/FG	440-2915	\$39.95

### CompoSite Finishing Kit (Shofu)

More shapes for a complete finishing system - 12 burs in each kit

Now in BurButler the universal gripping autoclavable silicone bur block

Use Composite Finishing kit for: composites

PN 0302B • Composite Finishing Metal Kit FG	440-2916	\$39.95
PN 0302 • Composite Finishing Plastic Kit FG	440-2917	\$39.95



### CompoSite Polishing Kit (Shofu)

More shapes for a complete finishing system - 14 burs in each kit

Now in BurButler the universal gripping autoclavable silicone bur block

Use Composite Finishing kit for: composites

PN 0310B • Composite Finishing Metal Kit CA	440-2913	\$39.95
PN 0310 • Composite Finishing Plastic Kit CA	440-2919	\$39.95



## Abrasives

kits

### CADMaster CAD/CAM™ (Shofu)

FG Kit 440-2900 \$114.95  
Dura-green stones (FL2 & 2 x RD2),  
Dura-white stones (FL2 & 2 x RD2),  
Ceramaster (2 x mini points),  
Ceramaster Coarse, Diamond Robot  
Point & T&F Hybrid Point (7406)



RA Kit 440-2902 \$146.95  
Dura-green stones (FL2, RD2, TC4 & WH6),  
Dura-white stones (FL2, RD2 & TC4),  
Ceramaster (bullet & WH6), Ceramaster  
Coarse (bullet & WH6)

HP Kit 440-2904 \$105.95  
Dura-green stones (FL2 & 2 x RD2), Dura-white  
stones (FL2 & 2 x RD2), Ceramaster (bullet,  
knife, cup & mini-point), Ceramaster Coarse  
(bullet, knife, cup & mini-point)

## Composite Polishing

interproximal

CPR Stain Remover (ICare)  
Perfect for composites, veneers  
and restorations, non-pumice  
100g cups/jar

Mint	440-5800	\$63.95
------	----------	---------





## Accessories

FG-RA

RA Mandrels (Microdont)

RA ID	Code# Size	ISO	Ø	Length mm	Description Details
	10502.001	611 - [323]	Ø23	13,0	FG Converted Adapted 13
	10502.002	603 - [303]	Ø60	22,0	Discs and Saws Discs & Saws
	10502.003	611 - [327]	Ø50	31,0	Polisher and Felt Polisher & Felt
	10502.004	603 - [303]	Ø30	19,0	Discs and Saws - Mini Discs & Saws - Mini
	10502.005	615 - [313]	Ø35	25,0	Discs - Pop-On Discs - Pop-On
	10502.006	611 - [327]	Ø45	23,0	Special Discs - Snap Special Discs - Snap
	10502.007	611 - [327]	Ø45	45,0	Special Discs - Snap Long Special Discs - Snap Long
	10502.009	612 - [327]	Ø50	23,0	Occlusal Polisher Occlusal Polisher

Disc Saw, 1/Pack	442-1100	\$0.95
Disc Saw, 1/Pack	442-1101	\$0.95
Pop-on Disc, 3/Pack	442-1102	\$14.95
Occlusal Polisher, 1/Pack	442-1103	\$0.95
Polisher Felt, 1/Pack	442-1104	\$0.95
Snap-on Disc, 1/Pack	442-1105	\$0.95
Snap-on Disc, 3/Pack	442-1106	\$5.95
FG Converted, 1/Pack	442-1107	\$4.95

For the latest in Composite  
Finishing & Polishing

It's Your Choice.

## Speedy

Economical, fast and easy snap-on  
silicone point system that produces  
a satin-like polish with one point

AVAILABLE IN 3 SNAP-ON SHAPES  
CONVENIENT STARTER KIT AVAILABLE  
REFILLS AVAILABLE IN PACKS OF 24

18 <sup>PIECES</sup>  
+  
MANDREL



## Poli

Utilizes a 3-grit sequence  
that works great on any composite.  
The combination produces  
superior esthetic results!

REFILLS AVAILABLE IN PACKS OF 24

18 <sup>PIECES</sup>  
+  
MANDREL



the Latest in Snap-On  
Silicone Systems



For the latest in Composite  
Finishing & Polishing

It's Your Choice.

## POLIDONT DIAMOND DISC SYSTEM

MICRODONT  
USA

The latest in Composite Disk Systems

UTILIZES  
DIAMOND PARTICLES  
TO ACHIEVE THE ULTIMATE RESULTS

THIN MYLAR DISCS  
PLUS MINI DISCS TO REACH  
DIFFICULT AREAS

3 EASY STEPS  
ONLY 3 COLOR-CODED GRITS  
TO ACHIEVE BEST RESULTS

CONVENIENT  
STARTER KIT AVAILABLE

REFILLS AVAILABLE  
IN PACKS OF 24



28 <sup>PIECES</sup>  
+  
MANDREL