



Dryz[®] blu

Gingival Hemostatic Retraction Paste

NEW!

Create a Dry Working Environment Every Time

- Bright aqua 'blu' color allows easy detection of where material is placed.
- Provides fast, effective hemostasis due to its aluminum chloride content.
- Stiff consistency for more displacement of sulcular tissue.
- Delivered in single dose capsules for more precise placement and asepsis.



NEW! Dryz Blu Intro Kit
Includes 30 Unit Dose Capsules and Dispensing Gun

NEW! Dryz Blu Unit Dose
Includes 30 Unit Dose Capsules



SHARP[™]

PREMIUM HYDROPHILIC VINYL POLYSILOXANE IMPRESSION MATERIAL



- Exceptional hydrophilicity: equal to or better than other leading brands.
- Outstanding dimensional accuracy for dead-on stone pours.
- Superb detail replication for well-fitting restorations every time.
- Excellent tear strength for easy removal without marginal tearing.



Color/ Viscosity	Shore-A Hardness
● Wash	52
● Medium	55
● Heavy	68

Includes two 50ml split cartridges (100ml ttl.) and mixing tips (Wash also includes intraoral tips.)



Follow Us!

parkell
PARKELL.COM

Retraction Cords

Gingi-Pak

Gingi-Plain® Soft-Twist
Purple Bottle, Non-impregnated

	108"	
	\$12.95	
#1 Thin	825-2711	
#2 Medium	825-2721	
#3 Large	825-2731	



Gingi-Pak® Soft-Twist
Green Bottle, w/Epinephrine

108" length	\$12.95
#1 Thin	825-2751
#2 Medium	825-2761
#3 Large	825-2771



Gingi-Plain Z-Twist
Purple Bottle, Non-impregnated

108" length	\$12.95
#00 X-Thin	825-2811
#1 Thin	825-2821
#2 Medium	825-2831
#3 Large	825-2841



Gingi-Pak Z-Twist
Green Bottle, w/Epinephrine

108" length	\$12.95
#00 X-Thin	825-2811
#1 Thin	825-2861
#2 Medium	825-2871
#3 Large	825-2881



Gingi-Aid Z-Twist
Blue Bottle, w/Aluminum Sulfate

108" length	\$12.95
#00 X-Thin	825-2895
#1 Thin	825-2896
#2 Medium	825-2897
#3 Large	825-2898



Dux/Kerr

Aluminium potassium sulfate
yellow cap

10 ft		\$15.95
#0a	825-8704	
#1a	825-8705	
#2a	825-8706	
#3a	825-8707	



Non-impregnated
white cap

10 ft		\$15.95
#0n	825-8708	
#1n	825-8709	
#2n	825-8710	
#3n	825-8711	



Unibraid
Unit-Dose GingiBraid, pre-cut &
individually packed in 2" lengths
50/pkg

Epinephrine/Aluminum, Blue			
#0e	825-8720	#2e	825-8722
#1e	825-8721	#3e	825-8723
Aluminium potassium sulfate, Yellow			
#0a	825-8724	#2a	825-8726
#1a	825-8725		



Gingiknit
Knitted Yarn, Bottle of 6 ft

Aluminum potassium sulfate Yellow cap			
#00a	825-8741	#1a	825-8743
#0a	825-8742	#2a	825-8744
		#3a	825-8745



Non-Impregnated white cap			
#00n	825-8747	#1n	825-8749
#0n	825-8748	#2n	825-8750
		#3n	825-8751

Fax your order to: 800.931.0601

RETRACTION MATERIALS

Retraction Cords

Pascal

Knit Trax
Plain knitted retraction cord
72"/ bottle

#00 826-6200
#0 826-6201
#1 826-6202
#2 826-6203



Racord®
W/Epinephrine Hydrochloride
72"/ bottle

#7 826-6207
#8 826-6208
#9 826-6209
#10 826-6210



Racord® II
W/Epinephrine Hydrochloride &
Zinc Phenol, 72"/ bottle

#7 826-6217
#8 826-6218
#9 826-6219
#10 826-6220



Premier

Knit Pak (Premier)
Non-impregnated, 100"/ bottle

#000 825-8815
#00 825-8816
#0 825-8817
#1 825-8818
#2 825-8819
#3 825-8820



Vista Dental

Knitted Stat-Pak
Color coded cords for easy identification,
Plain

100" length \$13.95



#000 Black 825-8900
#00 Blue 825-8901
#0 Purple 825-8902
#1 Light Blue 825-8903
#2 Yellow 825-8904

Stat-Pak + E
Braided retraction cord with
epinephrine,



72" length \$13.95

#00 Dark Purple 825-8910
#0 Light Purple 825-8911
#1 Light Blue 825-8912

Compression Caps

Parkell

Compression Caps

The Perfect Assist for the Perfect Impression.

Made of pure, naturally-absorbent cotton that has been specially-starched to maintain its shape during hard biting, Parkell's Compression Caps provide safe and effective gingival retraction and fluid absorption during deep, subgingival impressions.



Non-Anatomic Assortment #S196A 827-1100 \$39.95
150 total caps (3 different sizes, 50 of each size)

Anatomic Assortment #S195A 827-1101 \$42.95
120 total caps (3 different sizes, 40 of each size)

RETRACTION MATERIALS

Retraction Paste

Syringe

(Parkell)

Dryz® Blu Retraction Paste

Dryz® Blu features a contrasting blue color making the material more visible against tooth structure.



Syringe Kit 835-7016 \$69.95
7 Syringes + 15 tips

Dryz® Retraction Paste

Dryz® is terrific for use alone or in conjunction with retraction cord and compression caps.



Syringe Kit 836-2100 \$67.95
7 Syringes + 15 tips

Value Pack 836-2110 \$204.95
25 Syringes + 55 tips

Dryz® Luer Lock Applicator Tips 836-2115 \$31.95
40 tips

Application Tips

Sol-u-Flo™ (Vista)

For use with low viscosity solutions
19 ga, Opaque Black

20/cup #312202	830-0010	\$12.95
100/cup #312200	830-0011	\$43.95
500/cup #312205	830-0012	\$159.95



Stat-Flo™ (Vista)

For use with hemostatic agents
19 ga, Yellow

20-Pack #312102	830-0020	\$12.95
100-Pack #312100	830-0022	\$51.95



Retraction Paste

Capsules

Dryz® Blu Retraction Paste (Parkell)

Dryz® Blu features a contrasting blue color making the material more visible against tooth structure.



Intro Kit 835-7011 \$124.95
30 unit dose capsules (0.25 g/capsule) + Dispensing Gun

Standard Kit 835-7010 \$94.95
30 unit dose capsules (0.25 g/capsule)

Dryz® Retraction Paste

Dryz® is terrific for use alone or in conjunction with retraction cord and compression caps.



Intro Kit 835-7000 \$88.95
30 unit dose capsules (0.25 g/capsule) + Dispensing Gun

Hemostatic Solution

bottle

(Dharma)

Frenna Ac

25% aluminum chloride, retraction material, epinephrine free

1 oz	835-3005	\$12.95
------	----------	---------



Hemox-A (Keystone)

25% Aluminium Chloride

1oz	835-3000	\$12.95
-----	----------	---------



(Premier)

Hemodent

Buffered Aluminum Chloride, Epinephrine free

10cc	835-6210	\$14.95
20cc	835-6220	\$28.95
40cc	835-6240	\$51.95



RETRACTION MATERIALS

Hemostatic Solution

bottle

Hemoban (Sultan)
Aluminum Chloride, Epinephrine free
30ml 835-7500 \$43.95



Gingi-Aid (Gingi Pak)
25% buffered aluminium chloride
15ml 835-3500 \$15.95



Stasis (Gingi Pak)
21% ferric sulfide hemostyptic
15ml 835-3520 \$15.95

Orostat (Gingi Pak)
8% dl epinephrine HCl
15ml 835-3510 \$16.95

Rastringent II (Pascal)
with 25% Aluminum Sulfate
15ml 835-5012 \$14.95
40ml 835-5013 \$28.95



Quick-Stat™ Free (Vista)
Clear Hemostatic Agent. Does not stain,
controls minor bleeding fast, Stat-Flo™
padded tip

Standard Kit #504600 835-8120 \$19.95
4 x 1.2ml syringes + 8 x Stat-Flo tips



Bulk Syringe Kit #504650 835-8122 \$37.95
30ml filled syringe, 20 stat-flo tips,
20 x 1.2ml empty syringes, docking
port

Bulk Syringe Refill #504660 835-8125 \$28.95
30ml filled syringe, syringe docking port

Gel-Cord (Pascal)
25% Aluminum Sulfate
12 x 1ml syringes + tips 835-6300 \$42.95



Stat-Gel (Pascal)
15% Ferric Sulfate
12 x 1ml syringe + tips 835-6310 \$42.95



Hemostatic Solution

syringe

Quick-Stat™ FS (Vista)
stops bleeding when completing
impressions, restorations, crown &
bridge procedures, 15.5% Ferric Sulfate
Standard Kit #502804 835-8100 \$19.95
4 x 1.2ml syringes + 8 x Stat-Flo tips



Bulk Syringe Kit #503550 835-8110 \$37.95
30ml syringe + 20 x Stat-Flo tips +
20 x 1.2cc syringes

Bulk Syringe Refill #503555
30ml filled syringe, syringe docking port
835-8115 \$29.95



Gel Cord Clear (Pascal)
25% Aluminum Sulfate
12 x 1ml syringe + tips 835-6315 \$42.95



Racellet Hemostatic (Pascal)
Racemic Epinephrine HCl Pellets
12 x 1ml syringe + tips \$29.95



#2 Large 250/Pkg 835-6320
#3 Small 450/Pkg 835-6321

