






















Rubber Dam Clamps




stainless steel

(FCI)		#30	#31	#51	#138	#2	
Molar Clamps							
Regular		840-1156	840-1157	840-1158	840-1159	840-1160	
Each	\$4.99						
		W#3	W#4	#24	#25	#26	#28
							
		840-1161	840-1162	840-1163	840-1164	840-1165	840-1166
		#26A	#28A	#30A	#31A	#51A	
							
		840-1174	840-1175	840-1176	840-1177	840-1178	

(FCI)		#200	#201	#202	#203	#204	#205
Molar Clamps							
Winged		840-1182	840-1183	840-1184	840-1185	840-1186	840-1187
Each	\$4.99						

(FCI)		#27	#29	#8	#9	#10
Bicuspid Clamps						
Regular		840-1172	840-1173	840-1179	840-1180	840-1181
Each	\$4.99					

(FCI)		#206	#207	#208	#209	#22
Bicuspid Clamps						
Winged		840-1167	840-1168	840-1169	840-1170	840-1171
Each	\$4.99					

(FCI)		#210	#211	#212
Labial Clamps				
Winged		840-1188	840-1189	840-1190
Each	\$4.99			



RUBBER DAM

Rubber Dam Clamps

stainless steel

(Integra/Miltex)

Miltex rubber dam clamps are stainless steel and fully autoclavable.

Each \$10.49

Molars

Each \$10.49



Lower Molars



Upper Molars



Pre Molars

Each \$10.49




Pre Molars

Each \$10.49



Anterior

Each \$10.49



RUBBER DAM

Rubber Dam Clamps

stainless steel

Ivory Rubber Dam Clamp (Kulzer)




Complete Kit 840-0100 \$384.95
Contains punch, forceps, rubber dam, frame, 8 clamp starter kit & template


Starter Kit 840-0101 \$92.95
8 clamps - 0, 2A, 8A, 14, 2,7,9,14A

Individual Winged Clamps (Kulzer)


Anterior (Labial) \$16.95
Stainless Steel



Bicuspid \$16.95
Stainless Steel



Molar \$16.95
Stainless Steel



Rubber Dams

latex

Dental Dam (Blossom) \$6.95
Green Mint 5" x 5", 52 sheets 840-0001
6" x 6", 36 sheets 840-0002
Blue 5" x 5", 52 sheets 840-0003
6" x 6", 36 sheets 840-0004




Nic Tone Dental Dam (MDC)
Latex, low in powder & protein
5" x 5", 52 sheets
6" x 6", 36 sheets

	Thin	Medium	Heavy
Green 5x5	840-0010	840-0020	840-0030
Green 6x6	840-0011	840-0021	840-0031
Blue 5x5	840-0012	840-0022	840-0032
Blue 6x6	840-0013	840-0023	840-0033



IVORY® Rubber Dam

Made from natural latex, thoroughly washed in the manufacturing stage to reduce surface proteins. Ivory dam has extremely high tear strength, reducing waste from tearing during application. The dam is thin and pliable, allowing easy placement in especially difficult areas. The bright blue or green tint increases visibility by its contrast.



Available in 6x6 or 5x5 (medium weight).

IVORY® Punch

This precision instrument is made from high quality stainless steel. The unique floating plunger strikes all edges of the cutting disc evenly, guaranteeing a perfect, tear-resistant hole every time. The cutting disc rotates through a range of 6 hole sizes, providing the correct opening for the smallest anterior tooth to the largest molar.

IVORY® Forceps

The Ivory forceps are made from high grade stainless steel, which will assure many years of use. The forceps will reach around the bow of any clamp without tilting, which is especially important when reversing a clamp. The straight head design allows easy gripping and placement of any clamp size.

IVORY® Rubber Dam Kits

The Ivory Clamp Starter Kit offers a selection of the most widely used clamps. These popular clamps provide a variety of general use clamps that fulfill most needs. The Rubber Dam Complete Kit is a self-contained unit providing all the necessary equipment for the Rubber Dam Technique. The Complete Kit is especially useful for the practitioner just beginning their practice or for one who is implementing the technique for the first time.



RUBBER DAM

Rubber Dam

Latex

Fiesta Dam (Coltene)
3-color assortment pink, purple & blue, fruit scented
5" x 5", 52/pkg

840-0132 \$18.95

6" x 6", 36/pkg

840-0133 \$18.95



Regular Dam - Green (Coltene)

	Thin	Medium	Heavy
	\$14.95	\$15.95	\$17.95
5" x 5"	840-0130	840-0132	840-0134
6" x 6"	840-0131	840-0133	840-0135



Unflavored Blue

	Thin	Medium	Heavy
5" x 5"	840-0050	840-0055	840-0060
6" x 6"	840-0051	840-0056	840-0061

Mint Green

	Thin	Medium	Heavy
5" x 5"	840-0065	840-0070	840-0080
6" x 6"	840-0066	840-0075	840-0085

Smart Dam (Cranberry)
Latex, Medium

Spearmint, Aquagreen	\$11.95
5 x 5, 52/box	840-0144
6 x 6, 36/box	840-0145



Insti-Dams (Zirc)
Powder-free latex, 4 1/4" diameter, 20/pack
Natural
Dispenser

Natural	840-0150	\$49.95
Dispenser	840-0151	\$23.95



Dental Dam (Crosstex)
Latex Powder Free
\$11.95

Unflavored Blue

	Thin	Medium	Heavy
5" x 5"	840-0050	840-0055	840-0060
6" x 6"	840-0051	840-0056	840-0061

Mint Green

	Thin	Medium	Heavy
5" x 5"	840-0065	840-0070	840-0080
6" x 6"	840-0066	840-0075	840-0085



Rubber Dams

latex-free

Synthetic Rubber Dam (MDC)
Latex free, violet
15/box

5" x 5"	840-0040
6" x 6"	840-0041



Smart Dam (Cranberry)
15/box
Non-Latex, Medium

Spearmint Scented Purple Color
6" x 6" 840-0141



Dental Dam (Crosstex)
Latex free

6" x 6", 36 sheets \$19.95



Hygienic Dental Dam (Coltene)
Latex free

5" x 5", 52 sheets
6" x 6", 36 sheets



Price are subject to change.

ALL PICTURES SHOWN ARE FOR ILLUSTRATION PURPOSE ONLY. ACTUAL PRODUCT MAY VARY DUE TO PRODUCT ENHANCEMENT.



RUBBER DAM

Insti-Dam (Zirc)
Latex-free, 4 1/4"

Blue 840-0155 \$54.95
Dispenser 840-0156 \$23.95
holds 35 Insti-dams



Framed Non-Latex Flexi
Dam (Coltene)
10x10.5 cm diameter,
20/pack
Natural

840-0140 \$66.95



Punch Forceps

(First Choice) Each \$12.95



Ainsworth
17cm 550-0599




Ivory
17cm 550-0600

(Integra/Miltex)



Ainsworth 840-3000



IV Type 840-3010

Rubber Dam Instruments

Clamp Forceps

(FCI) Each \$17.95



Ivory Lightweight
17 cm 550-0596



Brewer
17.5cm 550-0597

Stoke
17.5cm 550-0598



(Integra/Miltex)

Lightweight Forceps

840-3052



Brewer 840-3050



University of
Washington Pattern 840-3051

Rubber Dam Frame

Plastic 4" x 4"

(Pulpdent) 840-0200 \$10.95



(Integra/Miltex)

Frame #76-25 840-0220 \$25.95

Stainless Steel
fully autoclavable

(First Choice) \$6.95

4" x 4" 840-0209
5" x 5" 840-0210
6" x 6" 840-0211



(Integra/Miltex)

Ready Dam Frame with pins 840-0221

